



### SPECIFICATION SHEET

**Product Name** : Xylenes  
**Product Number** : 990.019  
**Brand** : ISOLAB  
**CAS Number** : 1330-20-7  
**Formula** :  $C_6H_4(CH_3)_2$   
**Formula Weight** : 106.17 g/mol

TEST NAME	SPECIFICATION	UNIT
Purity (Xylene isomers and ethylbenzol, GC)	99 - 100 %	%
Ethyl benzene (GC)	≤ 25 %	%
Non-volatile matter	≤ 0.002 %	%
Infrared spectrum	passes test	

Quality Management  
Isolab Laborgerate GmbH

