



SPECIFICATION SHEET

Product Name : Zinc sulfate heptahydrate, Extra pure
Product Number : 996.032
Brand : ISOLAB
CAS Number : 7446-20-0
Formula : $\text{ZnSO}_4 \cdot 7 \text{H}_2\text{O}$
Formula Weight : 287.54 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay	99.0 – 108.7 %	%
Appearance	Colourless, transparent prisms.	
Solubility	Very soluble in water. Freely soluble in glycerin, insoluble in ethanol.	
Identification	Responds to the tests for zinc and sulphate	
Acidity	passes test	
As (Arsenic)	$\leq 0.0014 \%$	%
Pb (Lead)	$\leq 0.002 \%$	%
Alkali and alkaline earths	$\leq 0.9 \%$	%

Quality Management
Isolab Laborgerate GmbH

