



## SPECIFICATION SHEET

**Product Name** : Iodine sublimated, for analysis ACS,ISO,Reag. Ph Eur  
**Product Number** : 935.016  
**Brand** : ISOLAB  
**CAS Number** : 7553-56-2  
**Formula** : I<sub>2</sub>  
**Formula Weight** : 253.8 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay	≥ 99,8 %	%
Identity	passes test	
Bromide and Chlorine (as Cl)	≤ 0,005 %	%
Residue	≤ 0,01 %	%

Quality Management  
Isolab Laborgerate GmbH

This document has been produced electronically and is valid without a signature