



SPECIFICATION SHEET

Product Name : Ethanol 96 % , Extra pure
Product Number : 920.052
Brand : ISOLAB
CAS Number : 64-17-5
Formula : C₂H₅OH
Formula Weight : 46.07 g/mol

TEST NAME	SPECIFICATION	UNIT
Purity (GC)	≥ 95.0 %	mol/l
Density (d 20°C/20°C)	0.805 - 0.812 g/cm ³	g/cm ³
Colour(Pt-Co)	≤ 10	
Identity (IR)	passes test	
Appearance	clear colorless liquid	
Acidity or alkalinity	passes test	
UV absorption	passes test	
Volatile impurities (GC) (Methanol)	≤ 0.15 %	%
Water	≤ 5 %	%

Quality Management
Isolab Laborgerate GmbH

