



PRODUCT SPECIFICATION

Product Name : Acetone $\geq 99.5\%$, Pure
Product Number : 901.022
Brand : ISOLAB
CAS Number : 67-64-1
Formula : C_3H_6O
Formula Weight : 58.08 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (GC)	$\geq 99.5\%$	%
Acidity (as Acetic Acid)	$\leq 0.002\%$	%
Density (d 20°C/20°C)	0.790 - 0.794	
Color	≤ 5	Pt-co
Related substances (GC), Methanol (Impurity A)	$\leq 0.05\%$	%
Related substances (GC), 2-Propanol (Impurity B)	$\leq 0.05\%$	%
Related substances (GC), Benzene (Impurity C)	$\leq 0.0002\%$	%
Related substances (GC) (Other impurity)	$\leq 0.05\%$	%
Water	$\leq 0.4\%$	%

Quality Management
Isolab Laborgerate GmbH

