



### SPECIFICATION SHEET

**Product Name** : Agarose, Bioreagent, for molecular biology, low EEO  
**Product Number** : 902.012  
**Brand** : ISOLAB  
**CAS Number** : 9012-36-6  
**Formula** : -  
**Formula Weight** : -

TEST NAME	SPECIFICATION	UNIT
Moisture	≤ 10 %	%
Sulphated ash	≤ 0.75 %	%
Sulphate	≤ 0.2 %	%
Gel strength (1.5% Solution)	> 800 g/cm <sup>2</sup>	g/cm <sup>2</sup>
Gelling temperature	35 - 40 °C	°C
EEO	0.09 - 0.12	

Quality Management  
Isolab Laborgerate GmbH

