



SPECIFICATION SHEET

Product Name : Agarose, Bioreagent, for molecular biology, low EEO
Product Number : 902.012
Brand : ISOLAB
CAS Number : 9012-36-6
EC Number : 232-731-8
Formula : -
Formula Weight : -

TEST NAME	SPECIFICATION	UNIT
Moisture	≤ 10 %	%
Sulphated ash	≤ 0.75 %	%
Sulphate	≤ 0.2 %	%
Gel strength (1.5% Solution)	> 800 g/cm ²	g/cm ²
Gelling temperature	35 - 40 °C	°C
EEO	0.09 - 0.12	

Quality Management
Isolab Laborgerate GmbH

