



SPECIFICATION SHEET

Product Name : Aniline Blue (Spirit soluble), (C.I. 42775) for microscopy
Product Number : 903.D02
Brand : ISOLAB
CAS Number : 8004-91-9
Formula : $C_{32}H_{25}N_3O_9S_3Na_2$
Formula Weight : 737.74 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	Reddish-purple to black powder	
Absorption maxima (in Methanol)	591-595 nm	nm

Quality Management
Isolab Laborgerate GmbH

