



SPECIFICATION SHEET

Product Name : Bismarck brown Y (G), (C.I. 21000) for microscopy
Product Number : 906.D02
Brand : ISOLAB
CAS Number : 10114-58-6
Formula : $C_{18}H_{18}N_8 \cdot 2HCl$
Formula Weight : 419.31 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	Fine black powder	
Absorption maxima	About 457 nm	nm

Quality Management
Isolab Laborgerate GmbH

