



SPECIFICATION SHEET

Product Name : Buffer Solution pH 3.00 ± 0.02 at 20°C
(Potassium Hydrogen Phthalate, Hydrochloric Acid)
Product Number : 908.B03
Brand : ISOLAB
CAS Number : -
Formula : -
Formula Weight : -

| TEST NAME | SPECIFICATION | UNIT |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|------|
| pH | 2.98 - 3.02 at 20°C | - |
| Measured with a combination glass electrode after multiple point calibration with reference materials. It is certified traceable to the following National Institute of Standards and Technology (USA), SRM 189c Potassium Tetraoxalate and SRM 185i Potassium Hydrogen Phthalate. | | |

Quality Management
Isolab Laborgerate GmbH

