



SPECIFICATION SHEET

Product Name : Copper(II) sulfate pentahydrate, Extra pure
Product Number : 911.022
Brand : ISOLAB
CAS Number : 7758-99-8
Formula : $\text{CuO}_4\text{S} \cdot 5 \text{H}_2\text{O}$
Formula Weight : 249.68 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay	$\geq 98.0 \%$	%
Fe (Iron)	$\leq 0.01 \%$	%
Free acid (H_2SO_4)	$\leq 0.1 \%$	%
Contents of insoluble residue	$\leq 0.2 \%$	%

Quality Management
Isolab Laborgerate GmbH

