

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Mixture
Product name : Formaldehyde solution

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

No additional information available

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

ISOLAB Laborgeräte GmbH
Am Dillhof 2 - 63863 Eschau / GERMANY
Tel: + 49 93 74 / 978 55-0
Fax: +49 93 74 / 978 55-29
prodsafe@isolab.de

1.4. Emergency telephone number

Country	Organisation/Company	Address	Emergency number	Comment
Germany	Giftnotruf der Charité CBF, Haus VIII (Wirtschaftsgebäude), UG	Hindenburgdamm 30 12203 Berlin	+49 30 19240	

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Acute Tox. 3 (Oral) H301
Acute Tox. 3 (Dermal) H311
Acute Tox. 3 (Inhalation:dust,mist) H331
Skin Corr. 1B H314
Skin Sens. 1 H317
Muta. 2 H341
Carc. 1B H350
STOT SE 1 H370
STOT SE 3 H335

Full text of hazard classes and H-statements : see section 16

Adverse physicochemical, human health and environmental effects

May cause cancer. Suspected of causing genetic defects. May cause damage to organs. Toxic in contact with skin. Toxic if inhaled. Toxic if swallowed. May cause respiratory irritation. Causes severe skin burns and eye damage. May cause an allergic skin reaction.

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) :



GHS05

GHS06

GHS08

Signal word (CLP) :

Danger

Hazardous ingredients :

formaldehyde ...%; methanol

Hazard statements (CLP) :

H301+H311+H331 - Toxic if swallowed, in contact with skin or if inhaled
H314 - Causes severe skin burns and eye damage.
H317 - May cause an allergic skin reaction
H335 - May cause respiratory irritation.
H341 - Suspected of causing genetic defects.
H350 - May cause cancer.
H370 - Causes damage to organs.

Precautionary statements (CLP) :

P201 - Obtain special instructions before use.

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



P280 - Wear protective gloves/ protective clothing/ eye protection/ face protection.
P301 + P330 + P331 - IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P302 + P352 - IF ON SKIN: Wash with plenty of soap and water.
P304 + P340 - IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P308 + P310 - IF exposed or concerned: immediately call a POISON CENTER or doctor/ physician.

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
formaldehyde ... % (Note B)(Note D)	(CAS-No.) 50-00-0 (EC-No.) 200-001-8 (EC Index-No.) 605-001-00-5 (REACH-no) 01-2119488953-20	37 - 50	Carc. 1B, H350 Muta. 2, H341 Acute Tox. 3 (Oral), H301 Acute Tox. 3 (Dermal), H311 Acute Tox. 3 (Inhalation), H331 Skin Corr. 1B, H314 Skin Sens. 1, H317
methanol	(CAS-No.) 67-56-1 (EC-No.) 200-659-6 (EC Index-No.) 603-001-00-X	2 - 7	Flam. Liq. 2, H225 Acute Tox. 3 (Inhalation), H331 Acute Tox. 3 (Dermal), H311 Acute Tox. 3 (Oral), H301 STOT SE 1, H370
formic acid ... % (Note B)	(CAS-No.) 64-18-6 (EC-No.) 200-579-1 (EC Index-No.) 607-001-00-0	< 10	Skin Corr. 1A, H314

Specific concentration limits:

Name	Product identifier	Specific concentration limits
formaldehyde ... %	(CAS-No.) 50-00-0 (EC-No.) 200-001-8 (EC Index-No.) 605-001-00-5 (REACH-no) 01-2119488953-20	(C >= 0.2) Skin Sens. 1, H317 (C >= 5) STOT SE 3, H335 (5 =<C < 25) Eye Irrit. 2, H319 (5 =<C < 25) Skin Irrit. 2, H315 (C >= 25) Skin Corr. 1B, H314
methanol	(CAS-No.) 67-56-1 (EC-No.) 200-659-6 (EC Index-No.) 603-001-00-X	(3 =<C < 10) STOT SE 2, H371 (C >= 10) STOT SE 1, H370
formic acid ... %	(CAS-No.) 64-18-6 (EC-No.) 200-579-1 (EC Index-No.) 607-001-00-0	(2 =<C < 10) Eye Irrit. 2, H319 (2 =<C < 10) Skin Irrit. 2, H315 (10 =<C < 90) Skin Corr. 1B, H314 (C >= 90) Skin Corr. 1A, H314

Note B : Some substances (acids, bases, etc.) are placed on the market in aqueous solutions at various concentrations and, therefore, these solutions require different classification and labelling since the hazards vary at different concentrations. In Part 3 entries with Note B have a general designation of the following type: 'nitric acid ... %'. In this case the supplier must state the percentage concentration of the solution on the label. Unless otherwise stated, it is assumed that the percentage concentration is calculated on a weight/weight basis.

Note D : Certain substances which are susceptible to spontaneous polymerisation or decomposition are generally placed on the market in a stabilised form. It is in this form that they are listed in Part 3. However, such substances are sometimes placed on the market in a non-stabilised form. In this case, the supplier must state on the label the name of the substance followed by the words 'non-stabilised'.

Full text of H-statements: see section 16

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general : Call a physician immediately.
First-aid measures after inhalation : Remove person to fresh air and keep comfortable for breathing. Call a doctor.
First-aid measures after skin contact : Take off immediately all contaminated clothing. Rinse skin with water/shower. Call a physician immediately.

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



First-aid measures after eye contact : Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Call a physician immediately.

First-aid measures after ingestion : Rinse mouth. Call a physician immediately. Do not induce vomiting.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects after inhalation : May cause respiratory irritation.

Symptoms/effects after skin contact : Burns. May cause an allergic skin reaction.

Symptoms/effects after eye contact : Serious damage to eyes.

Symptoms/effects after ingestion : Burns.

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Water spray. Dry powder. Foam. Carbon dioxide.

5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in case of fire : Toxic fumes may be released.

5.3. Advice for firefighters

Protection during firefighting : Do not attempt to take action without suitable protective equipment. Self-contained breathing apparatus. Complete protective clothing.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Emergency procedures : Only qualified personnel equipped with suitable protective equipment may intervene. Do not breathe dust/fume/gas/mist/vapours/spray.

6.1.2. For emergency responders

Protective equipment : Do not attempt to take action without suitable protective equipment. For further information refer to section 8: "Exposure controls/personal protection".

6.2. Environmental precautions

Avoid release to the environment. Notify authorities if product enters sewers or public waters.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Take up liquid spill into absorbent material. Notify authorities if product enters sewers or public waters.

Other information : Dispose of materials or solid residues at an authorized site.

6.4. Reference to other sections

For further information refer to section 13.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Ensure good ventilation of the work station. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Take all necessary technical measures to avoid or minimize the release of the product on the workplace. Limit quantities of product at the minimum necessary for handling and limit the number of exposed workers. Provide local exhaust or general room ventilation. Wear personal protective equipment. Floors, walls and other surfaces in the hazard area must be cleaned regularly. Do not breathe dust/fume/gas/mist/vapours/spray. Do not get in eyes, on skin, or on clothing.

Hygiene measures : Separate working clothes from town clothes. Launder separately. Wash contaminated clothing before reuse. Contaminated work clothing should not be allowed out of the workplace. Do not eat, drink or smoke when using this product. Always wash hands after handling the product.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store locked up. Store in a well-ventilated place. Keep container tightly closed. Keep cool.

7.3. Specific end use(s)

No additional information available

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



formaldehyde ...% (50-00-0)		
EU	Local name	Formaldehyde
EU	IOELV TWA (ppm)	0.2 ppm
EU	IOELV STEL (ppm)	0.4 ppm
EU	Notes	skin sensitizer. SCOEL Recommendations (2008/Ongoing)
Austria	Local name	Formaldehyd
Austria	MAK (mg/m ³)	0.6 mg/m ³
Austria	MAK (ppm)	0.5 ppm
Austria	MAK Short time value (mg/m ³)	0.6 mg/m ³
Austria	MAK Short time value (ppm)	0.5 ppm
Austria	Remark (AT)	H,Sh
Belgium	Local name	Aldéhyde formique # Formaldehyde
Belgium	Short time value (mg/m ³)	0.38 mg/m ³
Belgium	Short time value (ppm)	0.3 ppm
Belgium	Remark (BE)	M: La mention M indique que lors d'une exposition supérieure à la valeur limite, des irritations apparaissent ou un danger d'intoxication aiguë existe. Le procédé de travail doit être conçu de telle façon que l'exposition ne dépasse jamais la valeur limite. Lors des mesurages, la période d'échantillonnage doit être aussi courte que possible afin de pouvoir effectuer des mesurages fiables. Le résultat des mesurages est calculé en fonction de la période d'échantillonnage.# De vermelding M duidt aan dat bij de blootstelling boven de grenswaarde irritatie optreedt of er gevaar bestaat voor acute vergiftiging. Het werkproces moet zo zijn ontworpen dat de blootstelling de grenswaarde nooit overschrijdt. Bij een controle geldt dat de bemonsterde periode zo kort mogelijk moet zijn om een betrouwbare meting te kunnen verrichten. het meetresultaat wordt dan gerelateerd aan de beschouwde periode.
Bulgaria	Local name	Формалдехид
Bulgaria	OEL TWA (mg/m ³)	1 mg/m ³
Bulgaria	OEL STEL (mg/m ³)	2 mg/m ³
Croatia	Local name	Formaldehyd
Croatia	GVI (granična vrijednost izloženosti) (mg/m ³)	2.5 mg/m ³
Croatia	GVI (granična vrijednost izloženosti) (ppm)	2 ppm
Croatia	KGVI (kratkotrajna granična vrijednost izloženosti) (mg/m ³)	2.5 mg/m ³
Croatia	KGVI (kratkotrajna granična vrijednost izloženosti) (ppm)	2 ppm
Croatia	Naznake (HR)	T (otrovno); Karc. kat. 3 (tvari koje izazivaju zabrinutost zbog mogućeg karcinogenog djelovanja na ljude)
Czech Republic	Local name	Formaldehyd
Czech Republic	Expoziční limity (PEL) (mg/m ³)	0.5 mg/m ³
Czech Republic	Expoziční limity (PEL) (ppm)	0.41 ppm
Czech Republic	Expoziční limity (NPK-P) (mg/m ³)	1 mg/m ³
Czech Republic	Expoziční limity (NPK-P) (ppm)	0.81 ppm
Czech Republic	Remark (CZ)	D, S
Denmark	Local name	Formaldehyd (Formalin)
Denmark	Grænseværdie (langvarig) (mg/m ³)	0.4 mg/m ³
Denmark	Grænseværdie (langvarig) (ppm)	0.3 ppm
Denmark	Anmærkninger (DK)	L (markerer, at grænseværdien er en loftværdi, som ikke på noget tidspunkt må overskrides); K (betyder, at stoffet anses for at kunne være kræftfremkaldende)
Finland	Local name	Formaldehydi

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



formaldehyde ...% (50-00-0)		
Finland	HTP-arvo (8h) (mg/m ³)	0.37 mg/m ³
Finland	HTP-arvo (8h) (ppm)	0.3 ppm
Finland	HTP-arvo (15 min)	1.2 mg/m ³
Finland	HTP-arvo (15 min) (ppm)	1 ppm
Finland	Huomautus (FI)	kattoarvo
France	Local name	Aldéhyde formique
France	VME (ppm)	0.5 ppm
France	VLE (ppm)	1 ppm
France	Note (FR)	Valeurs recommandées/admises; substance classée cancérigène de catégorie 2
Germany	Local name	Formaldehyd
Germany	TRGS 900 Occupational exposure limit value (mg/m ³)	0.37 mg/m ³
Germany	TRGS 900 Occupational exposure limit value (ppm)	0.3 ppm
Germany	Remark (TRGS 900)	AGS,Sh,Y,X
Greece	OEL TWA (mg/m ³)	2.5 mg/m ³
Greece	OEL TWA (ppm)	2 ppm
Greece	OEL STEL (mg/m ³)	2.5 mg/m ³
Greece	OEL STEL (ppm)	2 ppm
Hungary	Local name	FORMALDEHID
Hungary	AK-érték	0.6 mg/m ³
Hungary	CK-érték	0.6 mg/m ³
Hungary	Megjegyzések (HU)	b, m, sz; VI.
Ireland	Local name	Formaldehyde
Ireland	OEL (8 hours ref) (ppm)	0.2 ppm
Ireland	OEL (15 min ref) (ppm)	0.4 ppm
Ireland	Notes (IE)	Carc.1B, Sens.
Latvia	Local name	Formaldehyds (metanāls)
Latvia	OEL TWA (mg/m ³)	0.5 mg/m ³
Lithuania	Local name	Formaldehidas
Lithuania	IPRV (mg/m ³)	0.6 mg/m ³
Lithuania	IPRV (ppm)	0.5 ppm
Lithuania	NRV (mg/m ³)	1.2 mg/m ³
Lithuania	NRV (ppm)	1 ppm
Lithuania	Remark (LT)	Ū (ūmus poveikis); J (jautrinantis poveikis); K (kancerogeninis poveikis)
Netherlands	Local name	Formaldehyde
Netherlands	Grenswaarde TGG 8H (mg/m ³)	0.15 mg/m ³
Netherlands	Grenswaarde TGG 15MIN (mg/m ³)	0.5 mg/m ³
Poland	Local name	Formaldehyd
Poland	NDS (mg/m ³)	0.5 mg/m ³
Poland	NDSch (mg/m ³)	1 mg/m ³
Portugal	Local name	Formaldeído
Portugal	OEL - Ceilings (mg/m ³)	0.3 mg/m ³
Portugal	OEL - Ceilings (ppm)	0.3 ppm
Romania	Local name	Formaldehida
Romania	OEL TWA (mg/m ³)	1.2 mg/m ³
Romania	OEL TWA (ppm)	1 ppm
Romania	OEL STEL (mg/m ³)	3 mg/m ³
Romania	OEL STEL (ppm)	2 ppm
Slovakia	Local name	Formaldehyd (metanál)

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



formaldehyde ...% (50-00-0)		
Slovakia	NPHV (priemerná) (mg/m ³)	0.37 mg/m ³
Slovakia	NPHV (priemerná) (ppm)	0.3 ppm
Slovakia	OEL STEL (mg/m ³)	0.74 mg/m ³
Slovakia	OEL STEL (ppm)	0.6 ppm
Slovakia	Upozornenie (SK)	S - znamená, že faktor môže spôsobiť senzibilizáciu
Slovenia	Local name	formaldehyd
Slovenia	OEL TWA (mg/m ³)	0.62 mg/m ³
Slovenia	OEL TWA (ppm)	0.5 ppm
Slovenia	OEL STEL (mg/m ³)	0.62 mg/m ³
Slovenia	OEL STEL (ppm)	0.5 ppm
Slovenia	KTV factor SL	1
Spain	Local name	Formaldehído
Spain	VLA-EC (mg/m ³)	0.37 mg/m ³
Spain	VLA-EC (ppm)	0.3 ppm
Spain	Notes	Sen (Sensibilizante. Véase Apartado 6), s (Esta sustancia tiene prohibida total o parcialmente su comercialización y uso como fitosanitario y/o como biocida. Para una información detallada acerca de las prohibiciones consúltese: Base de datos de productos biocidas: http://www.msssi.gob.es/ciudadanos/productos.do?tip o=plaguicidas Base de datos de productos fitosanitarios http://www.magrama.gob.es/agricultura/pags/fitos/regi stro/fichas/pdf/Lista_sa.pdf).
Sweden	Local name	Formaldehyd
Sweden	nivågränsvärde (NVG) (mg/m ³)	0.37 mg/m ³
Sweden	nivågränsvärde (NVG) (ppm)	0.3 ppm
Sweden	kortidsvärde (KTV) (mg/m ³)	0.74 mg/m ³
Sweden	kortidsvärde (KTV) (ppm)	0.6 ppm
Sweden	Anmärkning (SE)	C (Ämnet är cancerframkallande Risk för cancer finns även vid annan exponering än via inandning. För vissa cancerframkallande ämnen som inte har gränsvärden gäller förbud eller tillståndskrav enligt föreskrifterna om kemiska arbetsmiljörisiker); H (Ämnet kan lätt upptas genom huden Det föreskrivna gränsvärdet bedöms ge tillräckligt skydd endast under förutsättning att huden är skyddad mot exponering för ämnet ifråga); M (Medicinska kontroller kan krävas för hantering av ämnet. Se vidare föreskrifterna om medicinska kontroller i arbetslivet. För vissa ämnen ska arbetsgivaren erbjuda läkarundersökning och för andra ämnen gäller krav på periodisk läkarundersökning och tjänstbarhetsbedömning. Se föreskrifterna om kemiska arbetsmiljörisiker); S (Ämnet är sensibiliserande Sensibiliserande ämnen kan ge allergi eller annan överkänslighet. Överkänslighetsbesvären drabbar främst huden eller andningsorganen. Överkänslighet innebär att man reagerar vid kontakt med ämnen som normalt inte ger besvär. Allergi är en undergrupp av överkänslighet som orsakas av reaktioner i kroppens immunsystem. Särskilt låga gränsvärden har fastställts för ämnen med mer uttalat luftvägssensibiliserande egenskaper. Några ämnen med starkt sensibiliserande egenskaper får endast hanteras efter tillstånd från Arbetsmiljöverket, se föreskrifterna om kemiska arbetsmiljörisiker. Dessa ämnen har inga gränsvärden men i vissa fall riktvärden)
United Kingdom	Local name	Formaldehyde
United Kingdom	WEL TWA (mg/m ³)	2.5 mg/m ³

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



formaldehyde ...% (50-00-0)		
United Kingdom	WEL TWA (ppm)	2 ppm
United Kingdom	WEL STEL (mg/m ³)	2.5 mg/m ³
United Kingdom	WEL STEL (ppm)	2 ppm
Russian Federation	Local name	Формальдегид+
Russian Federation	OEL Ceiling (mg/m ³)	0.5 mg/m ³
Russian Federation	Remark (RU)	2 класс опасности - высокоопасное; п (пары и/или газы); О (вещества с остронаправленным механизмом действия, требующие автоматического контроля за их содержанием в воздухе); А (вещества, способные вызывать аллергические заболевания в производственных условиях); + (соединения, при работе с которыми требуется специальная защита кожи и глаз; символ проставлен вслед за наименованием вещества)
Norway	Local name	Formaldehyd
Norway	Grenseverdier (AN) (mg/m ³)	0.6 mg/m ³ 1.2 mg/m ³
Norway	Grenseverdier (AN) (ppm)	0.5 ppm 1 ppm
Norway	Grenseverdier (Takverdi) (mg/m ³)	1.2 mg/m ³
Norway	Grenseverdier (Takverdi) (ppm)	1 ppm
Norway	Merknader (NO)	A (Kjemikalier som skal betraktes som at de fremkaller allergi eller annen overfølsomhet i øynene eller luftveier, eller som skal betraktes som at de fremkaller allergi ved hudkontakt); K (Kjemikalier som skal betraktes som kreftfremkallende)
Switzerland	Local name	Formaldehyd
Switzerland	MAK (mg/m ³)	0.37 mg/m ³
Switzerland	MAK (ppm)	0.3 ppm
Switzerland	KZGW (mg/m ³)	0.74 mg/m ³
Switzerland	KZGW (ppm)	0.6 ppm
Switzerland	Remark (CH)	S C1 _B # SS _C - Auge - DFG, HSE, NIOSH, OSHA, *Kein erhöhtes Krebsrisiko bei Einhalten des MAK-Wertes [§] . 1.3.2.3
Australia	Local name	Formaldehyde
Australia	TWA (mg/m ³)	1.2 mg/m ³
Australia	TWA (ppm)	1 ppm
Australia	Remark (AU)	Carcinogenicity Category 2 – Suspected human carcinogen. The classification of a chemical into this category is on the basis of evidence from human and animal studies, where the evidence is not sufficiently convincing to place the chemical into Category 1 or from limited evidence of carcinogenicity in human or animal studies; Sen - Respiratory and/or Skin Sensitiser.
USA - ACGIH	Local name	Formaldehyde
USA - ACGIH	ACGIH TWA (ppm)	0.1 ppm
USA - ACGIH	ACGIH STEL (ppm)	0.3 ppm
USA - ACGIH	Remark (ACGIH)	URT & eye irr; URT cancer; DSEN; RSEN; A1 (Confirmed Human Carcinogen: The agent is carcinogenic to humans based on the weight of evidence from epidemiologic studies)
methanol (67-56-1)		
EU	Local name	Methanol
EU	IOELV TWA (mg/m ³)	260 mg/m ³
EU	IOELV TWA (ppm)	200 ppm
EU	Notes	skin
Austria	Local name	Methanol
Austria	MAK (mg/m ³)	260 mg/m ³
Austria	MAK (ppm)	200 ppm

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



methanol (67-56-1)		
Austria	MAK Short time value (mg/m ³)	1040 mg/m ³
Austria	MAK Short time value (ppm)	800 ppm
Austria	Remark (AT)	H
Belgium	Local name	Alcool méthylique # Methanol
Belgium	Limit value (mg/m ³)	266 mg/m ³
Belgium	Limit value (ppm)	200 ppm
Belgium	Short time value (mg/m ³)	333 mg/m ³
Belgium	Short time value (ppm)	250 ppm
Belgium	Remark (BE)	D: La mention D signifie que la résorption de l'agent, via la peau, les muqueuses ou les yeux, constitue une partie importante de l'exposition totale. Cette résorption peut se faire tant par contact direct que par présence de l'agent dans l'air. # De vermelding D betekent dat de opname van het agens via de huid, de slijmvliezen of de ogen een belangrijk deel van de totale blootstelling vormt. Deze opname kan het gevolg zijn van zowel direct contact als zijn aanwezigheid in de lucht.
Bulgaria	Local name	Метилов алкохол
Bulgaria	OEL TWA (mg/m ³)	260 mg/m ³
Bulgaria	OEL TWA (ppm)	200 ppm
Bulgaria	Notes	Кожа (възможна е значителна резорбция чрез кожата); • (Химични агенти, за които са определени гранични стойности във въздуха на работната среда за Европейската общност)
Croatia	Local name	Metanol
Croatia	GVI (granična vrijednost izloženosti) (mg/m ³)	260 mg/m ³
Croatia	GVI (granična vrijednost izloženosti) (ppm)	200 ppm
Croatia	Naznake (HR)	K (Skin): (naznaka da tvar može štetno djelovati kroz kožu); EU** (naznaka da se radi o tvarima za koje su utvrđene indikativne granične vrijednosti izloženosti prema Direktivi 2006/15/ EC (druga lista)); F (lako zapaljivo); T (otrovno)
Czech Republic	Local name	Methanol
Czech Republic	Expoziční limity (PEL) (mg/m ³)	250 mg/m ³
Czech Republic	Expoziční limity (PEL) (ppm)	189 ppm
Czech Republic	Expoziční limity (NPK-P) (mg/m ³)	1000 mg/m ³
Czech Republic	Expoziční limity (NPK-P) (ppm)	750 ppm
Czech Republic	Remark (CZ)	D
Denmark	Local name	Methanol (Methylalkohol)
Denmark	Grænseværdie (langvarig) (mg/m ³)	260 mg/m ³
Denmark	Grænseværdie (langvarig) (ppm)	200 ppm
Denmark	Anmærkninger (DK)	E (betyder, at stoffet har en EF-grænseværdi); H (betyder, at stoffet kan optages gennem huden)
Estonia	Local name	Metanool (metüülalkohol)
Estonia	OEL TWA (mg/m ³)	250 mg/m ³
Estonia	OEL TWA (ppm)	200 ppm
Estonia	OEL STEL (mg/m ³)	350 mg/m ³
Finland	Local name	Metanoli
Finland	HTP-arvo (8h) (mg/m ³)	270 mg/m ³
Finland	HTP-arvo (8h) (ppm)	200 ppm
Finland	HTP-arvo (15 min)	330 mg/m ³
Finland	HTP-arvo (15 min) (ppm)	250 ppm
Finland	Huomautus (FI)	iho
France	Local name	Alcool méthylique (méthanol)

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



methanol (67-56-1)		
France	VME (mg/m ³)	260 mg/m ³
France	VME (ppm)	200 ppm
France	VLE (mg/m ³)	1300 mg/m ³
France	VLE (ppm)	1000 ppm
France	Note (FR)	VME réglementaires contraignantes; la VLE n'est pas réglementaire et provient d'une circulaire du ministère chargé du travail; risque de pénétration percutanée
Germany	Local name	Methanol
Germany	TRGS 900 Occupational exposure limit value (mg/m ³)	270 mg/m ³
Germany	TRGS 900 Occupational exposure limit value (ppm)	200 ppm
Germany	Remark (TRGS 900)	DFG,EU,H,Y
Greece	OEL TWA (mg/m ³)	260 mg/m ³
Greece	OEL TWA (ppm)	200 ppm
Greece	OEL STEL (mg/m ³)	325 mg/m ³
Greece	OEL STEL (ppm)	250 ppm
Hungary	Local name	METIL-ALKOHOL
Hungary	AK-érték	260 mg/m ³
Hungary	Megjegyzések (HU)	b, i; II.1.
Ireland	Local name	Methanol
Ireland	OEL (8 hours ref) (mg/m ³)	260 mg/m ³
Ireland	OEL (8 hours ref) (ppm)	200 ppm
Ireland	Notes (IE)	Sk, IOELV
Italy	Local name	Metanolo
Italy	OEL TWA (mg/m ³)	260 mg/m ³
Italy	OEL TWA (ppm)	200 ppm
Latvia	Local name	Metanols (metilspirts, karbinols)
Latvia	OEL TWA (mg/m ³)	260 mg/m ³
Latvia	OEL TWA (ppm)	200 ppm
Lithuania	Local name	Metanolis (metilo alkoholis)
Lithuania	IPRV (mg/m ³)	260 mg/m ³
Lithuania	IPRV (ppm)	200 ppm
Lithuania	Remark (LT)	O (medžiaga į organizmą gali prasiskverbti pro nepažeistą odą)
Luxembourg	Local name	Méthanol
Luxembourg	OEL TWA (mg/m ³)	260 mg/m ³
Luxembourg	OEL TWA (ppm)	200 ppm
Malta	Local name	Methanol
Malta	OEL TWA (mg/m ³)	260 mg/m ³
Malta	OEL TWA (ppm)	200 ppm
Netherlands	Local name	Methanol
Netherlands	Grenswaarde TGG 8H (mg/m ³)	133 mg/m ³
Netherlands	Remark (MAC)	H (Huidopname) Stoffen die relatief gemakkelijk door de huid kunnen worden opgenomen, hetgeen een substantiële bijdrage kan betekenen aan de totale inwendige blootstelling, hebben in de lijst een Haanduiding. Bij deze stoffen moeten naast maatregelen tegen inademing ook adequate maatregelen ter voorkoming van huidcontact worden genomen.
Poland	Local name	Metanol (metylowy alkohol)
Poland	NDS (mg/m ³)	100 mg/m ³
Poland	NDSch (mg/m ³)	300 mg/m ³
Portugal	Local name	Metanol (Álcool metílico)

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



methanol (67-56-1)		
Portugal	OEL TWA (ppm)	200 ppm
Portugal	OEL STEL (ppm)	250 ppm
Romania	Local name	Alcool metilic
Romania	OEL TWA (mg/m ³)	260 mg/m ³
Romania	OEL TWA (ppm)	200 ppm
Romania	OEL STEL (ppm)	5 ppm
Slovakia	Local name	Metylalkohol (metanol)
Slovakia	NPHV (priemerná) (mg/m ³)	260 mg/m ³
Slovakia	NPHV (priemerná) (ppm)	200 ppm
Slovakia	Upozornenie (SK)	K - znamená, že faktor môže byť ľahko absorbovaný kožou
Slovenia	Local name	metanol (metilalkohol)
Slovenia	OEL TWA (mg/m ³)	260 mg/m ³
Slovenia	OEL TWA (ppm)	200 ppm
Spain	Local name	Metanol (Alcohol metílico)
Spain	VLA-ED (mg/m ³)	266 mg/m ³
Spain	VLA-ED (ppm)	200 ppm
Spain	Notes	Vía dérmica (Indica que, en las exposiciones a esta sustancia, la aportación por la vía cutánea puede resultar significativa para el contenido corporal total si no se adoptan medidas para prevenir la absorción. En estas situaciones, es aconsejable la utilización del control biológico para poder cuantificar la cantidad global absorbida del contaminante. Para más información véase el Apartado 5 de este documento), VLB® (Agente químico que tiene Valor Límite Biológico específico en este documento), VLI (Agente químico para el que la U.E. estableció en su día un valor límite indicativo. Todos estos agentes químicos figuran al menos en una de las directivas de valores límite indicativos publicadas hasta ahora (ver Anexo C. Bibliografía). Los estados miembros disponen de un tiempo fijado en dichas directivas para su transposición a los valores límites de cada país miembro. Una vez adoptados, estos valores tienen la misma validez que el resto de los valores adoptados por el país).
Sweden	Local name	Metanol
Sweden	nivågränsvärde (NVG) (mg/m ³)	250 mg/m ³
Sweden	nivågränsvärde (NVG) (ppm)	200 ppm
Sweden	kortidsvärde (KTV) (mg/m ³)	350 mg/m ³
Sweden	kortidsvärde (KTV) (ppm)	250 ppm
Sweden	Anmärkning (SE)	H (Ämnet kan lätt upptas genom huden Det föreskrivna gränsvärdet bedöms ge tillräckligt skydd endast under förutsättning att huden är skyddad mot exponering för ämnet ifråga); V (Vägledande korttidsgränsvärde ska användas som ett rekommenderat högsta värde som inte bör överskridas)
United Kingdom	Local name	Methanol
United Kingdom	WEL TWA (mg/m ³)	266 mg/m ³
United Kingdom	WEL TWA (ppm)	200 ppm
United Kingdom	WEL STEL (mg/m ³)	333 mg/m ³
United Kingdom	WEL STEL (ppm)	250 ppm
United Kingdom	Remark (WEL)	Sk (Can be absorbed through the skin. The assigned substances are those for which there are concerns that dermal absorption will lead to systemic toxicity)
Russian Federation	Local name	Метанол+
Russian Federation	OEL Ceiling (mg/m ³)	15 mg/m ³

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



methanol (67-56-1)		
Russian Federation	OEL TWA (mg/m ³)	5 mg/m ³
Russian Federation	Remark (RU)	3 класс опасности - опасное; п (пары и/или газы); + (соединения, при работе с которыми требуется специальная защита кожи и глаз; символ проставлен вслед за наименованием вещества)
Norway	Local name	Metanol
Norway	Grenseverdier (AN) (mg/m ³)	130 mg/m ³
Norway	Grenseverdier (AN) (ppm)	100 ppm
Norway	Merknader (NO)	H (Kjemikalier som kan tas opp gjennom huden); E (EU har en veiledende grenseverdi for stoffet)
Switzerland	Local name	Methanol
Switzerland	MAK (mg/m ³)	260 mg/m ³ 260 mg/m ³
Switzerland	MAK (ppm)	200 ppm 200 ppm
Switzerland	KZGW (mg/m ³)	1040 mg/m ³ 1040 mg/m ³
Switzerland	KZGW (ppm)	800 ppm 800 ppm
Switzerland	Remark (CH)	H B SS _c - ZNS, Sehen - INRS, NIOSH
Turkey	Local name	Metanol
Turkey	OEL TWA (mg/m ³)	260 mg/m ³
Turkey	OEL TWA (ppm)	200 ppm
Turkey	Comments	Deri
Australia	Local name	Methyl alcohol
Australia	TWA (mg/m ³)	262 mg/m ³ Synonym (Methanol)
Australia	TWA (ppm)	200 ppm Synonym (Methanol)
Australia	STEL (mg/m ³)	328 mg/m ³ Synonym (Methanol)
Australia	STEL (ppm)	250 ppm Synonym (Methanol)
Australia	Remark (AU)	Sk - Absorption through the skin may be a significant source of exposure.
USA - ACGIH	Local name	Methanol
USA - ACGIH	ACGIH TWA (ppm)	200 ppm
USA - ACGIH	ACGIH STEL (ppm)	250 ppm
USA - ACGIH	Remark (ACGIH)	Headache; eye dam; dizziness; nausea
USA - OSHA	Local name	Methyl alcohol
USA - OSHA	OSHA PEL (TWA) (mg/m ³)	260 mg/m ³
USA - OSHA	OSHA PEL (TWA) (ppm)	200 ppm
formic acid ... % (64-18-6)		
EU	Local name	Formic acid
EU	IOELV TWA (mg/m ³)	9 mg/m ³
EU	IOELV TWA (ppm)	5 ppm
Austria	Local name	Ameisensäure
Austria	MAK (mg/m ³)	9 mg/m ³
Austria	MAK (ppm)	5 ppm
Austria	MAK Short time value (mg/m ³)	9 mg/m ³
Austria	MAK Short time value (ppm)	5 ppm
Belgium	Local name	Acide formique # Mierenzuur
Belgium	Limit value (mg/m ³)	9.5 mg/m ³
Belgium	Limit value (ppm)	5 ppm
Belgium	Short time value (mg/m ³)	19 mg/m ³
Belgium	Short time value (ppm)	10 ppm
Bulgaria	Local name	Мравчена киселина
Bulgaria	OEL TWA (mg/m ³)	9 mg/m ³
Bulgaria	OEL TWA (ppm)	5 ppm

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



formic acid ... % (64-18-6)		
Bulgaria	Notes	• (Химични агенти, за които са определени гранични стойности във въздуха на работната среда за Европейската общност)
Croatia	Local name	Mravlja kiselina, > 90%
Croatia	GVI (granična vrijednost izloženosti) (mg/m ³)	9 mg/m ³
Croatia	GVI (granična vrijednost izloženosti) (ppm)	5 ppm
Croatia	Naznake (HR)	EU** (naznaka da se radi o tvarima za koje su utvrđene indikativne granične vrijednosti izloženosti prema Direktivi 2006/15/ EC (druga lista)); C (nagrizajuće)
Czech Republic	Local name	Kyselina mravení
Czech Republic	Expoziční limity (PEL) (mg/m ³)	9 mg/m ³
Czech Republic	Expoziční limity (PEL) (ppm)	4.78 ppm
Czech Republic	Expoziční limity (NPK-P) (mg/m ³)	18 mg/m ³
Czech Republic	Expoziční limity (NPK-P) (ppm)	9.56 ppm
Denmark	Local name	Myresyre (Methansyre)
Denmark	Grænseværdie (langvarig) (mg/m ³)	9 mg/m ³
Denmark	Grænseværdie (langvarig) (ppm)	5 ppm
Denmark	Anmærkninger (DK)	E (betyder, at stoffet har en EF-grænseværdi)
Estonia	Local name	Metaanhape (sipelghape)
Estonia	OEL TWA (mg/m ³)	9 mg/m ³
Estonia	OEL TWA (ppm)	5 ppm
Finland	Local name	Muurahaishappo
Finland	HTP-arvo (8h) (mg/m ³)	5 mg/m ³
Finland	HTP-arvo (8h) (ppm)	3 ppm
Finland	HTP-arvo (15 min)	19 mg/m ³
Finland	HTP-arvo (15 min) (ppm)	10 ppm
France	Local name	Acide formique
France	VME (mg/m ³)	9 mg/m ³
France	VME (ppm)	5 ppm
France	Note (FR)	Valeurs réglementaires indicatives
Germany	Local name	Ameisensäure
Germany	TRGS 900 Occupational exposure limit value (mg/m ³)	9.5 mg/m ³
Germany	TRGS 900 Occupational exposure limit value (ppm)	5 ppm
Germany	Remark (TRGS 900)	DFG,EU,Y
Greece	OEL TWA (mg/m ³)	9 mg/m ³
Greece	OEL TWA (ppm)	5 ppm
Hungary	Local name	HANGYASAV
Hungary	AK-érték	9 mg/m ³
Hungary	Megjegyzések (HU)	m; l.
Ireland	Local name	Formic acid
Ireland	OEL (8 hours ref) (mg/m ³)	9 mg/m ³
Ireland	OEL (8 hours ref) (ppm)	5 ppm
Ireland	Notes (IE)	IOELV
Italy	Local name	Acido formico
Italy	OEL TWA (mg/m ³)	9 mg/m ³
Italy	OEL TWA (ppm)	5 ppm
Latvia	Local name	Skudrskābe (metānskābe)
Latvia	OEL TWA (mg/m ³)	9 mg/m ³
Latvia	OEL TWA (ppm)	5 ppm

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



formic acid ... % (64-18-6)		
Lithuania	Local name	Skruzdzių rūgštis
Lithuania	IPRV (mg/m ³)	9 mg/m ³
Lithuania	IPRV (ppm)	5 ppm
Luxembourg	Local name	Acide formique
Luxembourg	OEL TWA (mg/m ³)	9 mg/m ³
Luxembourg	OEL TWA (ppm)	5 ppm
Malta	Local name	Formic acid
Malta	OEL TWA (mg/m ³)	9 mg/m ³
Malta	OEL TWA (ppm)	5 ppm
Netherlands	Local name	Mierenzuur
Netherlands	Grenswaarde TGG 15MIN (mg/m ³)	5 mg/m ³
Poland	Local name	Kwas mrówkowy
Poland	NDS (mg/m ³)	5 mg/m ³
Poland	NDSch (mg/m ³)	15 mg/m ³
Portugal	Local name	Ácido fórmico
Portugal	OEL TWA (ppm)	5 ppm
Portugal	OEL STEL (ppm)	10 ppm
Romania	Local name	Acid formic
Romania	OEL TWA (mg/m ³)	9 mg/m ³
Romania	OEL TWA (ppm)	5 ppm
Slovakia	Local name	Kyselina mravčia (kyselina metánová)
Slovakia	NPHV (priemerná) (mg/m ³)	9 mg/m ³
Slovakia	NPHV (priemerná) (ppm)	5 ppm
Slovenia	Local name	mravljična kislina
Slovenia	OEL TWA (mg/m ³)	9 mg/m ³
Slovenia	OEL TWA (ppm)	5 ppm
Spain	Local name	Ácido fórmico
Spain	VLA-ED (mg/m ³)	9 mg/m ³
Spain	VLA-ED (ppm)	5 ppm
Spain	Notes	VLI (Agente químico para el que la U.E. estableció en su día un valor límite indicativo. Todos estos agentes químicos figuran al menos en una de las directivas de valores límite indicativos publicadas hasta ahora (ver Anexo C. Bibliografía). Los estados miembros disponen de un tiempo fijado en dichas directivas para su transposición a los valores límites de cada país miembro. Una vez adoptados, estos valores tienen la misma validez que el resto de los valores adoptados por el país), s (Esta sustancia tiene prohibida total o parcialmente su comercialización y uso como fitosanitario y/o como biocida. Para una información detallada acerca de las prohibiciones consúltese: Base de datos de productos biocidas: http://www.msssi.gob.es/ciudadanos/productos.do?tipo=plaguicidas Base de datos de productos fitosanitarios: http://www.magrama.gob.es/agricultura/pags/fitos/registro/fichas/pdf/Lista_sa.pdf)
Sweden	Local name	Myrsyra
Sweden	nivågränsvärde (NVG) (mg/m ³)	5 mg/m ³
Sweden	nivågränsvärde (NVG) (ppm)	3 ppm
Sweden	kortidsvärde (KTV) (mg/m ³)	9 mg/m ³
Sweden	kortidsvärde (KTV) (ppm)	5 ppm

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



formic acid ... % (64-18-6)		
Sweden	Anmärkning (SE)	V (Vägledande korttidsgränsvärde ska användas som ett rekommenderat högsta värde som inte bör överskridas)
United Kingdom	Local name	Formic acid
United Kingdom	WEL TWA (mg/m ³)	9.6 mg/m ³
United Kingdom	WEL TWA (ppm)	5 ppm
Russian Federation	Local name	Метановая кислота+
Russian Federation	OEL Ceiling (mg/m ³)	1 mg/m ³
Russian Federation	Remark (RU)	2 класс опасности - высокоопасное; п (пары и/или газы); + (соединения, при работе с которыми требуется специальная защита кожи и глаз; символ проставлен вслед за наименованием вещества)
Norway	Local name	Maursyre
Norway	Grenseverdier (AN) (mg/m ³)	9 mg/m ³
Norway	Grenseverdier (AN) (ppm)	5 ppm
Norway	Merknader (NO)	E (EU har en veiledende grenseverdi for stoffet)
Switzerland	Local name	Ameisensäure
Switzerland	MAK (mg/m ³)	9.5 mg/m ³
Switzerland	MAK (ppm)	5 ppm
Switzerland	KZGW (mg/m ³)	19 mg/m ³
Switzerland	KZGW (ppm)	10 ppm
Switzerland	Remark (CH)	SS _C - Haut & Auge, OAW ^{KT AN} - NIOSH, OSHA
Turkey	Local name	Formik asit
Turkey	OEL TWA (mg/m ³)	9 mg/m ³
Turkey	OEL TWA (ppm)	5 ppm
Australia	Local name	Formic acid
Australia	TWA (mg/m ³)	9.4 mg/m ³
Australia	TWA (ppm)	5 ppm
Australia	STEL (mg/m ³)	19 mg/m ³
Australia	STEL (ppm)	10 ppm
USA - ACGIH	Local name	Formic acid
USA - ACGIH	ACGIH TWA (ppm)	5 ppm
USA - ACGIH	ACGIH STEL (ppm)	10 ppm
USA - ACGIH	Remark (ACGIH)	URT, eye, & skin irr
USA - OSHA	Local name	Formic acid
USA - OSHA	OSHA PEL (TWA) (mg/m ³)	9 mg/m ³
USA - OSHA	OSHA PEL (TWA) (ppm)	5 ppm

8.2. Exposure controls

Appropriate engineering controls:

Ensure good ventilation of the work station.

Hand protection:

Protective gloves

Eye protection:

Protective goggles

Skin and body protection:

Wear suitable protective clothing

Respiratory protection:

[In case of inadequate ventilation] wear respiratory protection.

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



Environmental exposure controls:

Avoid release to the environment.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state	: Liquid
Appearance	: Colourless liquid.
Colour	: Colourless.
Odour	: characteristic.
Odour threshold	: No data available
pH	: 2.8 - 4.0 (H ₂ O, 20 °C)
Relative evaporation rate (butylacetate=1)	: No data available
Melting point	: < -15 °C
Freezing point	: No data available
Boiling point	: 93 - 96 °C (1013 hPa)
Flash point	: 69 °C water solution.
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Flammability (solid, gas)	: Not applicable
Vapour pressure	: > 100 (kPa)
Relative vapour density at 20 °C	: No data available
Density	: 1.09 g/cm ³ (20 °C)
Solubility	: Miscible with water.
Log Pow	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive properties	: No data available
Oxidising properties	: No data available
Explosive limits	: 7 - 73 %(V) (Formaldehyde)

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No dangerous reactions known under normal conditions of use.

10.4. Conditions to avoid

None under recommended storage and handling conditions (see section 7).

10.5. Incompatible materials

No additional information available

10.6. Hazardous decomposition products

Under normal conditions of storage and use, hazardous decomposition products should not be produced.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Oral: Toxic if swallowed. Dermal: Toxic in contact with skin. Inhalation:dust,mist: Toxic if inhaled.

ATE CLP (oral)	175.4385964912 mg/kg bodyweight
ATE CLP (dermal)	526.3157894737 mg/kg bodyweight
ATE CLP (dust,mist)	0.8771929825 mg/l/4h

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



formaldehyde ...% (50-00-0)	
LD50 oral	100 mg/kg bodyweight
LD50 dermal	270 mg/kg bodyweight
LC50 inhalation rat (Dust/Mist - mg/l/4h)	497 mg/l

Skin corrosion/irritation	: Causes severe skin burns and eye damage.
Serious eye damage/irritation	: Serious eye damage, category 1, implicit
Respiratory or skin sensitisation	: May cause an allergic skin reaction.
Germ cell mutagenicity	: Suspected of causing genetic defects.
Carcinogenicity	: May cause cancer.
Reproductive toxicity	: Not classified
STOT-single exposure	: May cause damage to organs. May cause respiratory irritation.
STOT-repeated exposure	: Not classified
Aspiration hazard	: Not classified

SECTION 12: Ecological information

12.1. Toxicity

Ecology - general : Before neutralisation, the product may represent a danger to aquatic organisms.

formaldehyde ...% (50-00-0)	
LC50 fish 1	6.7 mg/l
EC50 other aquatic organisms 1	5.8 mg/l waterflea
EC50 other aquatic organisms 2	3.48 mg/l

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

formaldehyde ...% (50-00-0)	
Log Pow	0.779

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

No additional information available

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods : Dispose of contents/container in accordance with licensed collector's sorting instructions.

SECTION 14: Transport information

In accordance with ADR / RID / IMDG / IATA / ADN

ADR	IMDG	IATA	ADN	RID
14.1. UN number				
2209	2209	2209	2209	2209
14.2. UN proper shipping name				
FORMALDEHYDE SOLUTION (CONTAINS ; formaldehyde ... %)	FORMALDEHYDE SOLUTION (CONTAINS ; formaldehyde ... %)	Formaldehyde solution (CONTAINS ; formaldehyde ... %)	FORMALDEHYDE SOLUTION (CONTAINS ; formaldehyde ... %)	FORMALDEHYDE SOLUTION (CONTAINS ; formaldehyde ... %)
Transport document description				
UN 2209 FORMALDEHYDE SOLUTION (CONTAINS ; formaldehyde ... %), 8, III, (E)	UN 2209 FORMALDEHYDE SOLUTION (CONTAINS ; formaldehyde ... %), 8, III	UN 2209 Formaldehyde solution (CONTAINS ; formaldehyde ... %), 8, III	UN 2209 FORMALDEHYDE SOLUTION (CONTAINS ; formaldehyde ... %), 8, III	UN 2209 FORMALDEHYDE SOLUTION (CONTAINS ; formaldehyde ... %), 8, III

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2




ADR	IMDG	IATA	ADN	RID
14.3. Transport hazard class(es)				
8	8	8	8	8
14.4. Packing group				
III	III	III	III	III
14.5. Environmental hazards				
Dangerous for the environment : No	Dangerous for the environment : No Marine pollutant : No	Dangerous for the environment : No	Dangerous for the environment : No	Dangerous for the environment : No
No supplementary information available				

14.6. Special precautions for user

- Overland transport

Classification code (ADR)	: C9
Special provisions (ADR)	: 533
Limited quantities (ADR)	: 5I
Excepted quantities (ADR)	: E1
Packing instructions (ADR)	: P001, IBC03, LP01, R001
Mixed packing provisions (ADR)	: MP19
Portable tank and bulk container instructions (ADR)	: T4
Portable tank and bulk container special provisions (ADR)	: TP1
Tank code (ADR)	: L4BN
Vehicle for tank carriage	: AT
Transport category (ADR)	: 3
Special provisions for carriage - Packages (ADR)	: V12
Hazard identification number (Kemler No.)	: 80
Orange plates	:



Tunnel restriction code (ADR) : E

- Transport by sea

Limited quantities (IMDG)	: 5 L
Excepted quantities (IMDG)	: E1
Packing instructions (IMDG)	: P001, LP01
IBC packing instructions (IMDG)	: IBC03
Tank instructions (IMDG)	: T4
Tank special provisions (IMDG)	: TP1
EmS-No. (Fire)	: F-A
EmS-No. (Spillage)	: S-B
Stowage category (IMDG)	: A
Properties and observations (IMDG)	: Colourless, clear liquid, with a suffocating pungent odour. Usually stabilized with methyl alcohol. Miscible with water. Causes burns to skin, eyes and mucous membranes.
MFAG-No	: 132

- Air transport

PCA Excepted quantities (IATA)	: E1
PCA Limited quantities (IATA)	: Y841
PCA limited quantity max net quantity (IATA)	: 1L
PCA packing instructions (IATA)	: 852

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



PCA max net quantity (IATA) : 5L
CAO packing instructions (IATA) : 856
CAO max net quantity (IATA) : 60L
ERG code (IATA) : 8i

- Inland waterway transport

Classification code (ADN) : C9
Special provisions (ADN) : 533
Limited quantities (ADN) : 5 L
Excepted quantities (ADN) : E1
Carriage permitted (ADN) : T
Equipment required (ADN) : PP, EP
Number of blue cones/lights (ADN) : 0

- Rail transport

Classification code (RID) : C9
Special provisions (RID) : 533
Limited quantities (RID) : 5L
Excepted quantities (RID) : E1
Packing instructions (RID) : P001, IBC03, LP01, R001
Mixed packing provisions (RID) : MP19
Portable tank and bulk container instructions (RID) : T4
Portable tank and bulk container special provisions (RID) : TP1
Tank codes for RID tanks (RID) : L4BN
Transport category (RID) : 3
Special provisions for carriage – Packages (RID) : W12
Colis express (express parcels) (RID) : CE8
Hazard identification number (RID) : 80

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Contains no REACH substances with Annex XVII restrictions
Contains no substance on the REACH candidate list
Contains no REACH Annex XIV substances

15.1.2. National regulations

Germany

VwVwS Annex reference : Water hazard class (WGK) 2, significant hazard to waters (Classification according to VwVwS, Annex 4)

12th Ordinance Implementing the Federal Immission Control Act - 12.BImSchV : Is not subject of the 12. BImSchV (Hazardous Incident Ordinance)

Netherlands

SZW-lijst van kankerverwekkende stoffen : formaldehyde ...% is listed
SZW-lijst van mutagene stoffen : None of the components are listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding : None of the components are listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Vruchtbaarheid : None of the components are listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling : methanol is listed

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



Denmark

Class for fire hazard	: Class III-1
Store unit	: 50 liter
Classification remarks	: Flammable according to the Danish Ministry of Justice; Emergency management guidelines for the storage of flammable liquids must be followed
Recommendations Danish Regulation	: Young people below the age of 18 years are not allowed to use the product Pregnant/breastfeeding women working with the product must not be in direct contact with the product The requirements from the Danish Working Environment Authorities regarding work with carcinogens must be followed during use and disposal

15.2. Chemical safety assessment

No chemical safety assessment has been carried out

SECTION 16: Other information

Abbreviations and acronyms:

ADN	European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
ATE	Acute Toxicity Estimate
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
PBT	Persistent Bioaccumulative Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
RID	Regulations concerning the International Carriage of Dangerous Goods by Rail
SDS	Safety Data Sheet
vPvB	Very Persistent and Very Bioaccumulative

Data sources : REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

Full text of H- and EUH-statements:	
Acute Tox. 3 (Dermal)	Acute toxicity (dermal), Category 3
Acute Tox. 3 (Inhalation)	Acute toxicity (inhal.), Category 3
Acute Tox. 3 (Inhalation:dust,mist)	Acute toxicity (inhalation:dust,mist) Category 3
Acute Tox. 3 (Oral)	Acute toxicity (oral), Category 3
Carc. 1B	Carcinogenicity, Category 1B
Flam. Liq. 2	Flammable liquids, Category 2
Muta. 2	Germ cell mutagenicity, Category 2
Skin Corr. 1A	Skin corrosion/irritation, Category 1A
Skin Corr. 1B	Skin corrosion/irritation, Category 1B
Skin Sens. 1	Skin sensitisation, Category 1
STOT SE 1	Specific target organ toxicity — single exposure, Category 1
STOT SE 2	Specific target organ toxicity — Single exposure, Category 2
STOT SE 3	Specific target organ toxicity — Single exposure, Category 3, Respiratory tract irritation
H225	Highly flammable liquid and vapour.
H301	Toxic if swallowed.
H311	Toxic in contact with skin.
H314	Causes severe skin burns and eye damage
H317	May cause an allergic skin reaction
H331	Toxic if inhaled.
H335	May cause respiratory irritation.
H341	Suspected of causing genetic defects.
H350	May cause cancer.
H370	Causes damage to organs.

Formaldehyde solution

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 01/06/2022 Version: 0.0

Doc. No: SDS-923.012/2



SDS ISOLAB

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product