



### SPECIFICATION SHEET

Product Name : Formaldehyde solution about 37%, GR for analysis ACS, Reag. Ph Eur  
Product Number : 923.015  
Brand : ISOLAB  
CAS Number : -  
Formula : -  
Formula Weight : -

TEST NAME	SPECIFICATION	UNIT
Assay (Acidimetric)	36.5 - 38.0 %	%
Assay (Iodometric)	36.5 - 38.0 %	
Identity	passes test	
Colour	≤ 10 Hazen	
Appearance of solution	passes test	
Free acid (as HCOOH)	≤ 0.025 %	%
Density (20°C)	1.080 - 1.090 g/ml	g/ml
Chloride (Cl)	≤ 0.0001 %	%
Sulphate (SO <sub>4</sub> )	≤ 0.002 %	%
Heavy metals (as Pb)	≤ 0.0002 %	%
Fe (Iron)	≤ 0.0001 %	%
Methanol (GC)	9.0 - 11.0 %	%
Sulfated ash	≤ 0.002 %	%

Quality Management  
Isolab Laborgerate GmbH

