



SPECIFICATION SHEET

Product Name : Gallic acid monohydrate
Product Number : 925.011
Brand : ISOLAB
CAS Number : 5995-86-8
Formula : $C_7H_6O_5 \cdot H_2O$
Formula Weight : 188.13 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay	$\geq 99.0 \%$	%
Loss on drying	$\leq 10 \%$	%
Residue on ignition	$\leq 0.1 \%$	%
Sulphate	$\leq 0.02 \%$	%
Tannic acid	$\leq 1 \%$	%
Turbidity	$\leq 0.0010 \%$	%

Quality Management
Isolab Laborgerate GmbH

