



SPECIFICATION SHEET

Product Name : n-Heptane, Extra pure
Product Number : 930.012
Brand : ISOLAB
CAS Number : 142-82-5
Formula : C₇H₁₆
Formula Weight : 100.2 g/mol

TEST NAME	SPECIFICATION	UNIT
Purity (GC)	≥ 99.0 %	%
Identity (IR)	passes test	
Free acid (as CH ₃ COOH)	≤ 0.005 %	%
Density (d 20°C/ 4°C)	0.683 - 0.684	g/cm ³
Matter discoloured by H ₂ SO ₄	≤ 10 Hazen	Hazen
Evaporation residue	≤ 0.005 %	%

Quality Management
Isolab Laborgerate GmbH

