



### SPECIFICATION SHEET

**Product Name** : Hydrogen peroxide 35 %, Extra pure  
**Product Number** : 932.122  
**Brand** : ISOLAB  
**CAS Number** : -  
**Formula** : -  
**Formula Weight** : -

TEST NAME	SPECIFICATION	UNIT
Assay	34.5 - 36.0 %	%
Appearance	Clear colourless liquid	
Density (g/cm <sup>3</sup> ) 25 °C	1.128 - 1.132	g/cm3
Color (APHA)	< 10	
Conductivity	> 50 µS/sec	µS/sec
Acidity	≤ 0.03 %	%
Stability	≤ 0.3 ml/min	ml/min
Chloride (Cl <sup>-</sup> )	≤ 0.0001 %	%
Nitrate (NO <sub>3</sub> )	≤ 0.015 %	%
Phosphate (PO <sub>4</sub> )	≤ 0.005 %	%
Sulfate (SO <sub>4</sub> )	≤ 0.001 %	%
As (Arsenic)	≤ 0.0003 %	%
Cu (Copper)	≤ 0.0001 %	%
Fe (Iron)	≤ 0.00005 %	%
Mn (Manganese)	≤ 0.0001 %	%
Pb (Lead)	≤ 0.001 %	%
Sn (Tin)	≤ 0.001 %	%
Non volatile matter	≤ 0.015 %	%

Quality Management  
Isolab Laborgerate GmbH

