



### SPECIFICATION SHEET

**Product Name** : Iodine sublimated, for analysis ACS, ISO, Reag. Ph Eur  
**Product Number** : 935.016  
**Brand** : ISOLAB  
**CAS Number** : 7553-56-2  
**Formula** : I<sub>2</sub>  
**Formula Weight** : 253.8 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay	≥ 99.8 %	%
Identity	passes test	
Bromide and chlorine (as Cl)	≤ 0.005 %	%
Residue	≤ 0.01 %	%

Quality Management  
Isolab Laborgerate GmbH

