



SPECIFICATION SHEET

Product Name : Mercury(II) sulfate, for analysis ACS
Product Number : 947.036
Brand : ISOLAB
CAS Number : 7783-35-9
Formula : HgSO_4
Formula Weight : 296.65 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (complexometric)	$\geq 99.0 \%$	%
Chloride (Cl)	$\leq 0.003 \%$	%
Nitrate (NO_3)	$\leq 0.005 \%$	%
Ca (Calcium)	$\leq 0.005 \%$	%
Fe (Iron)	$\leq 0.005 \%$	%
Mercury (I) (as Hg)	$\leq 0.15 \%$	%
Mg (Magnesium)	$\leq 0.005 \%$	%
Na (Sodium)	$\leq 0.005 \%$	%
Residue after reduction	$\leq 0.02 \%$	%

Quality Management
Isolab Laborgerate GmbH

