



**SPECIFICATION SHEET**

**Product Name** : Methanol Dry, for Karl Fischer titration  
**Product Number** : 947.080  
**Brand** : ISOLAB  
**CAS Number** : 67-56-1  
**Formula** : CH<sub>4</sub>O  
**Formula Weight** : 32.04 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay	≥ 99.5 %	%
Appearance	A clear liquid	
Suitability (for KF-titration)	passes test	
Water (KF)	≤ 0.02 %	%

Quality Management  
Isolab Laborgerate GmbH

