



### SPECIFICATION SHEET

**Product Name** : Methylene chloride, Extra pure  
**Product Number** : 947.092  
**Brand** : ISOLAB  
**CAS Number** : 75-09-2  
**Formula** : CH<sub>2</sub>Cl<sub>2</sub>  
**Formula Weight** : 84.93 g/mol

TEST NAME	SPECIFICATION	UNIT
Purity	≥ 99.5 %	%
Appearance	Clear	
Density	1.31 - 1.33 g/cm <sup>3</sup>	g/cm <sup>3</sup>
Acidity	≤ 0.0005 meq/gr	meq/gr
Water (K.F)	≤ 0.1 %	%
Colour(Pt-Co)	≤ 10	

Quality Management  
Isolab Laborgerate GmbH

