



SPECIFICATION SHEET

Product Name : Methyl violet, (C.I. 42535) for microscopy
Product Number : 947.D03
Brand : ISOLAB
CAS Number : 8004-87-3
Formula : -
Formula Weight : -

TEST NAME	SPECIFICATION	UNIT
Appearance	Dark green crystalline powder	
Absorption maxima in ethanol	583 - 587 nm	nm
Absorptivity (1%, 1 cm, λ_{\max})	≥ 1600	
Ratio $\lambda_{\max} \pm 15$ nm	1.060 - 1.190	
Loss on drying (at 110 °C)	≤ 10 %	%

Quality Management
Isolab Laborgerate GmbH

