



### SPECIFICATION SHEET

**Product Name** : New fuchsin, (C.I. 42520) for microscopy  
**Product Number** : 951.D01  
**Brand** : ISOLAB  
**CAS Number** : 3248-91-7  
**Formula** :  $C_{22}H_{24}ClN_3$   
**Formula Weight** : 365.91 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	Dark green powder	
Assay	≥ 75 %	%
Absorption maxima (50% ethanol)	548 - 558 nm	nm

Quality Management  
Isolab Laborgerate GmbH

