



### SPECIFICATION SHEET

**Product Name** : Orcein, For microscopy  
**Product Number** : 957.D02  
**Brand** : ISOLAB  
**CAS Number** : 1400-62-0  
**Formula** : -  
**Formula Weight** : -

TEST NAME	SPECIFICATION	UNIT
Appearance	Dark brown to black-purplish powder	
Absorption in NaOH 0.01 mol/L	575 - 582 nm	nm
Specific absorptivity	≥ 500	
Loss on drying	≤ 7 %	%

Quality Management  
Isolab Laborgerate GmbH

