



SPECIFICATION SHEET

Product Name : Patent blue V (VF), (C.I.42045) for microscopy
Product Number : 958.D01
Brand : ISOLAB
CAS Number : 129-17-9
Formula : $C_{27}H_{31}N_2NaO_6S_2$
Formula Weight : 566.67 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	Dark Blue to Very Dark Blue and Dark Blue-Green to Very Dark Blue-Green	
Absorption maxima	637 - 639 nm	nm

Quality Management
Isolab Laborgerate GmbH

