



SPECIFICATION SHEET

Product Name : Phenol, GR for analysis ACS, Reag. Ph Eur.
Product Number : 959.045
Brand : ISOLAB
CAS Number : 108-95-2
Formula : C₆H₆O
Formula Weight : 94.11 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (bromometric)	99.0 - 100.5 %	%
Identity	passes test	
Appearance of solution (67 g/l, water)	clear and not more intensely coloured than reference solution B6	
Acidic substances	passes test	
Solidification temperature	≥ 39.5 °C	°C
Residue on evaporation	≤ 0.05 %	%
Water (according to Karl Fischer)	≤ 0.5 %	%

Quality Management
Isolab Laborgerate GmbH

