



SPECIFICATION SHEET

Product Name : PHENOL $\geq 99\%$, CRYSTAL,GR FOR ANALYSIS ACS,REAG. PH EUR
Product Number : 959.046
Brand : ISOLAB
CAS Number : 108-95-2
Formula : C_6H_6O
Formula Weight : 94.11 g/mol

TEST NAME	SPECIFICATION	UNIT
Purity (GC)	$\geq 99.0 \%$	%
Identity	passes test	
Appearance of solution (67 g/l, water)	clear and not more intensely soloured than reference solution B6	
Acidic substances	passes test	
Solidification temperature	$\geq 39.5 \text{ }^{\circ}\text{C}$	$^{\circ}\text{C}$
Residue on evaporation	$\leq 0.05 \%$	%
Water (according to Karl Fischer)	$\leq 0.5 \%$	%

Quality Management
Isolab Laborgerate GmbH

