



SPECIFICATION SHEET

Product Name : Potassium chromate, For analysis ACS, Reag. Ph Eur
Product Number : 960.046
Brand : ISOLAB
CAS Number : 7789-00-6
Formula : K_2CrO_4
Formula Weight : 194.19 g/mol

| TEST NAME | SPECIFICATION | UNIT |
|-----------------------------|----------------|------|
| Assay | 99.5 - 100.5 % | % |
| In water insoluble matter | ≤ 0.003 % | % |
| pH-value (5%; water) | 9.0 - 9.8 | |
| Chloride (Cl) | ≤ 0.001 % | % |
| Sulphate (SO ₄) | ≤ 0.01 % | % |
| Ca (Calcium) | ≤ 0.005 % | % |
| Na (Sodium) | ≤ 0.02 % | % |
| Pb (Lead) | ≤ 0.005 % | % |
| Loss on drying | ≤ 0.5 % | % |

Quality Management
Isolab Laborgerate GmbH

