



SPECIFICATION SHEET

Product Name : Potassium dichromate, For analysis ACS, ISO, Reag. Ph Eur
Product Number : 960.056
Brand : ISOLAB
CAS Number : 7778-50-9
Formula : $K_2Cr_2O_7$
Formula Weight : 294.19 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (Iodometric, calc. on dried substance)	≥ 99.9 %	%
Insoluble matter	≤ 0.003 %	%
Chloride (Cl)	≤ 0.001 %	%
Sulphate (SO ₄)	≤ 0.005 %	%
Ca (Calcium)	≤ 0.002 %	%
Cu (Copper)	≤ 0.001 %	%
Fe (Iron)	≤ 0.001 %	%
Na (Sodium)	≤ 0.02 %	%
Pb (Lead)	≤ 0.005 %	%
Loss on drying (105 °C)	≤ 0.05 %	%

Quality Management
Isolab Laborgerate GmbH

