



SPECIFICATION SHEET

Product Name : Polysorbate 80, For synthesis
Product Number : 960.161
Brand : ISOLAB
CAS Number : 9005-65-6
Formula : -
Formula Weight : -

| TEST NAME | SPECIFICATION | UNIT |
|-----------------------|---------------|-------------------|
| Density (d 20°C/ 4°C) | 1.060 - 1.090 | g/cm ³ |
| Saponification value | 45 - 55 | mg KOH/g |
| Hydroxyl value | 65 - 80 | mg KOH/g |
| Identity (IR) | passes test | |

Due to its specific melting range the product may be solid, liquid, a solidified melt or a supercooled melt.

Quality Management
Isolab Laborgerate GmbH

