



SPECIFICATION SHEET

Product Name : Potassium Permanganate, Extra pure
Product Number : 960.172
Brand : ISOLAB
CAS Number : 7722-64-7
Formula : KMnO_4
Formula Weight : 158.03 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (manganometric)	99.0 - 100.5 %	%
Identity	passes test	
Appearance of solution	passes test	
In water insoluble matter	≤ 0.2 %	%
Chloride (Cl)	≤ 0.01 %	%
Sulfate (SO_4)	≤ 0.03 %	%

Quality Management
Isolab Laborgerate GmbH

