



SPECIFICATION SHEET

Product Name : 1-Propanol, Extra pure
Product Number : 961.012
Brand : ISOLAB
CAS Number : 71-23-8
Formula : C₃H₈O
Formula Weight : 60.1 g/mol

TEST NAME	SPECIFICATION	UNIT
Purity (GC)	≥ 99.0 %	%
Identity (IR)	passes test	
Acidity (C ₂ H ₅ COOH)	≤ 0.002 %	%
Methanol (GC)	≤ 0.1 %	%
Readily carbonizable substances	passes test	
Water	≤ 0.20 %	%

Quality Management
Isolab Laborgeräte GmbH

