



PRODUCT SPECIFICATION

Product Name	:	2-Propanol $\geq 99.8\%$ Extra Pure
Product Number	:	961.022
Brand	:	ISOLAB
CAS Number	:	67-63-0
Formula	:	C_3H_8O
Formula Weight	:	60.1 g/mol

TEST NAME	SPECIFICATION	UNIT
Purity	≥ 99.8	%
Water	≤ 0.1	%
Acidity as acetic acid	≤ 0.001	%
Water Miscibility	Miscible	-
Nonvolatile Residue	≤ 0.001	g/100ml
Appearance	Clear	-
Color, Pt-Co	Max 5.0	Pt Co
Density 20°C	0.785 - 0.786	g/ml
Distillation, IBP	Min. 81.8	°C
Distillation, DP	Max. 82.8	°C

Quality Management
Isolab Laborgerate GmbH

