



SPECIFICATION SHEET

Product Name : Salicylic acid, Extra pure
Product Number : 967.012
Brand : ISOLAB
CAS Number : 69-72-7
Formula : $C_7H_6O_3$
Formula Weight : 138.12 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance	White, crystalline, powder	
Assay	≥ 99.5 %	%
Residue on ignition	≤ 0.1 %	%
Heavy metals	≤ 0.002 %	%
Chloride (Cl)	≤ 0.01 %	%
Sulfate	≤ 0.03 %	%
Loss on drying	≤ 0.5 %	%
Phenol	≤ 0.05 %	%

Quality Management
Isolab Laborgerate GmbH

