

Sodium acetate, trihydrate

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 11/12/2016

Doc No: SDS-969.016/1



SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Substance
Substance name : Sodium acetate, trihydrate
EC no : 204-823-8
CAS No : 6131-90-4
Type of product : Pure substance
Formula : C₂H₃NaO₂·3H₂O
Chemical structure :



Synonyms : acetic acid sodium salt, trihydrate / acetic acid, sodium salt, trihydrate / sodium acetate trihydrate
BIG no : 16464

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Use of the substance/mixture : No data available

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

ISOLAB Laborgeräte GmbH
Am Dillhof 2 - 63863 Eschau / GERMANY
Tel: + 49 93 74 / 978 55-0
Fax: +49 93 74 / 978 55-29
prodsafe@isolab.de

1.4. Emergency telephone number

Country	Organisation/Company	Address	Emergency number	Comment
Germany	Giftnotruf der Charité CBF, Haus VIII (Wirtschaftsgebäude), UG	Hindenburgdamm 30 12203 Berlin	+49 (0) 30 19240	

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not classified

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

EUH-statements : EUH210 - Safety data sheet available on request

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substances

Name	Product identifier	%
Sodium acetate, trihydrate	(CAS No) 6131-90-4 (EC no) 204-823-8	99 - 100

Full text of H-statements: see section 16

Sodium acetate, trihydrate

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 11/12/2016

Doc No: SDS-969.016/1



3.2. Mixtures

Not applicable

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general	: If you feel unwell, seek medical advice.
First-aid measures after inhalation	: Remove the victim into fresh air. Respiratory problems: consult a doctor/medical service.
First-aid measures after skin contact	: Rinse with water. Soap may be used. Take victim to a doctor if irritation persists.
First-aid measures after eye contact	: Rinse with water. Do not apply neutralizing agents. Take victim to an ophthalmologist if irritation persists.
First-aid measures after ingestion	: Rinse mouth with water. Immediately after ingestion: give lots of water to drink. Victim is fully conscious: immediately induce vomiting. Call Poison Information Centre (www.big.be/antigif.htm). Consult a doctor/medical service if you feel unwell.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/injuries after inhalation	: EXPOSURE TO HIGH CONCENTRATIONS: Coughing. Slight irritation.
Symptoms/injuries after skin contact	: Slight irritation.
Symptoms/injuries after eye contact	: Slight irritation.
Symptoms/injuries after ingestion	: AFTER ABSORPTION OF HIGH QUANTITIES: Nausea. Vomiting. Abdominal pain.
Chronic symptoms	: ON CONTINUOUS/REPEATED EXPOSURE/CONTACT: Red skin. Dry skin. Skin rash/inflammation. Runny nose. Irritation of the respiratory tract. Irritation of the nasal mucous membranes.

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media	: Water spray. Polyvalent foam. Alcohol-resistant foam. ABC powder. Carbon dioxide.
Unsuitable extinguishing media	: Container may slop over if solid jet (water/foam) is applied.

5.2. Special hazards arising from the substance or mixture

Fire hazard	: DIRECT FIRE HAZARD. Not easily combustible. In finely divided state: increased fire hazard. INDIRECT FIRE HAZARD. Heating increases the fire hazard. Reactions involving a fire hazard: see "Reactivity Hazard".
Explosion hazard	: DIRECT EXPLOSION HAZARD. Fine dust is explosive with air. INDIRECT EXPLOSION HAZARD. Dust cloud can be ignited by a spark.

5.3. Advice for firefighters

Precautionary measures fire	: Exposure to fire/heat: keep upwind. Exposure to fire/heat: consider evacuation. Exposure to fire/heat: have neighbourhood close doors and windows.
Firefighting instructions	: Cool tanks/drums with water spray/remove them into safety. Dilute toxic gases with water spray.
Protection during firefighting	: Heat/fire exposure: compressed air/oxygen apparatus.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Protective equipment	: Gloves. Protective clothing. Dust cloud production: compressed air/oxygen apparatus. Reactivity hazard: compressed air/oxygen apparatus. Reactivity hazard: gas-tight suit. See "Material-Handling" to select protective clothing.
Emergency procedures	: Mark the danger area. Prevent dust cloud formation, e.g. by wetting. No naked flames. Wash contaminated clothes. In case of hazardous reactions: keep upwind. In case of reactivity hazard: consider evacuation.
Measures in case of dust release	: In case of dust production: keep upwind. Dust production: have neighbourhood close doors and windows. Dust production: stop engines and no smoking. In case of dust production: no naked flames or sparks. Dust: spark-/explosionproof appliances/lighting equipment.

6.1.2. For emergency responders

No additional information available

6.2. Environmental precautions

No additional information available

Sodium acetate, trihydrate

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 11/12/2016

Doc No: SDS-969.016/1



6.3. Methods and material for containment and cleaning up

- For containment : Contain released substance, pump into suitable containers. Plug the leak, cut off the supply. Knock down/dilute dust cloud with water spray. Provide equipment/receptacles with earthing. Powdered form: no compressed air for pumping over spills. If reacting: dilute toxic gas/vapour with water spray. Take account of toxic/corrosive precipitation water.
- Methods for cleaning up : Prevent dust cloud formation. Scoop solid spill into closing containers. Powdered: do not use compressed air for pumping over spills. Contaminated surfaces: clean (treat) with an excess of water. Wash clothing and equipment after handling.

6.4. Reference to other sections

No additional information available

SECTION 7: Handling and storage

7.1. Precautions for safe handling

- Precautions for safe handling : Comply with the legal requirements. Remove contaminated clothing immediately. Clean contaminated clothing. Thoroughly clean/dry the installation before use. Powdered form: no compressed air for pumping over. Avoid raising dust. Take precautions against electrostatic charges. Keep away from naked flames/heat. Finely divided: spark- and explosionproof appliances. Finely divided: keep away from ignition sources/sparks. Use earthed equipment. Observe normal hygiene standards. Keep container tightly closed. Measure the concentration in the air regularly. Carry operations in the open/under local exhaust/ventilation or with respiratory protection.

7.2. Conditions for safe storage, including any incompatibilities

- Heat and ignition sources : KEEP SUBSTANCE AWAY FROM: heat sources. ignition sources.
- Information on mixed storage : KEEP SUBSTANCE AWAY FROM: oxidizing agents. (strong) acids. water/moisture.
- Storage area : Store in a dry area. Keep container in a well-ventilated place. Provide the tank with earthing. Meet the legal requirements.
- Special rules on packaging : SPECIAL REQUIREMENTS: closing. watertight. dry. clean. correctly labelled. meet the legal requirements. Secure fragile packagings in solid containers.
- Packaging materials : SUITABLE MATERIAL: No data available. MATERIAL TO AVOID: No data available.

7.3. Specific end use(s)

No additional information available

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Sodium acetate, trihydrate (6131-90-4)		
Belgium	Limit value (mg/m ³)	3 mg/m ³ (Particules non classifiées autrement (fraction alvéolaire); Belgium; Time-weighted average exposure limit 8 h; Particules non classifiées autrement (fraction inhalable); 10 mg/m ³ ; Belgium; Time-weighted average exposure limit 8 h)
France	VME (mg/m ³)	10 mg/m ³ (Poussières réputées sans effet spécifique; France; Time-weighted average exposure limit 8 h; VRC: Valeur réglementaire contraignante; Poussières réputées sans effet spécifique, fraction; 5 mg/m ³ ; France; Time-weighted average exposure limit 8 h; VRC: Valeur réglementaire contraignante)
United Kingdom	WEL TWA (mg/m ³)	4 mg/m ³ Inhalable dust; United Kingdom; Time-weighted average exposure limit 8 h; Workplace exposure limit (EH40/2005); Respirable dust; 10 mg/m ³ ; United Kingdom; Time-weighted average exposure limit 8 h; Workplace exposure limit (EH40/2005)

8.2. Exposure controls

Materials for protective clothing:

GIVE EXCELLENT RESISTANCE: No data available. GIVE GOOD RESISTANCE: butyl rubber. PVC. GIVE LESS RESISTANCE: No data available. GIVE POOR RESISTANCE: No data available

Hand protection:

Gloves

Eye protection:

Sodium acetate, trihydrate

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 11/12/2016

Doc No: SDS-969.016/1



Safety glasses. In case of dust production: protective goggles

Skin and body protection:

Protective clothing

Respiratory protection:

Dust production: dust mask with filter type P1

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state	: Solid
Appearance	: Crystalline solid. Powder. Grains. Solid.
Molecular mass	: 136.08 g/mol
Colour	: Colourless.
Odour	: Odourless.
Odour threshold	: No data available
pH	: 7.5 - 9.2 (50 g/l, H ₂ O, 20 °C)
pH solution	: 5 %
Relative evaporation rate (butylacetate=1)	: No data available
Melting point	: 58 °C
Freezing point	: No data available
Boiling point	: Not applicable
Flash point	: > 250 °C
Auto-ignition temperature	: > 600 °C
Decomposition temperature	: 58 °C
Flammability (solid, gas)	: No data available
Vapour pressure	: No data available
Relative vapour density at 20 °C	: No data available
Relative density	: 1.5
Density	: 1450 kg/m ³
Solubility	: Soluble in water. Soluble in ethanol. Soluble in ether. Water: > 61.3 g/100ml Ethanol: 5.2 g/100ml
Log Pow	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive properties	: No data available
Oxidising properties	: No data available
Explosive limits	: No data available

9.2. Other information

VOC content	: 0 %
Other properties	: Substance has basic reaction.

SECTION 10: Stability and reactivity

10.1. Reactivity

Decomposes on exposure to temperature rise: release of corrosive gases/vapours (acetic acid vapours). Upon combustion: CO and CO₂ are formed. Reacts violently with (strong) oxidizers. Reacts violently with (some) acids: release of corrosive gases/vapours (acetic acid vapours).

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No additional information available

10.4. Conditions to avoid

No additional information available

10.5. Incompatible materials

No additional information available

Sodium acetate, trihydrate

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 11/12/2016

Doc No: SDS-969.016/1



10.6. Hazardous decomposition products

No additional information available

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Not classified

Sodium acetate, trihydrate (6131-90-4)	
LD50 oral rat	3530 mg/kg
LD50 dermal rabbit	> 10000 mg/kg (Rabbit)

Skin corrosion/irritation : Not classified
pH: 7.5 - 9.0 (5 %)

Serious eye damage/irritation : Not classified
pH: 7.5 - 9.0 (5 %)

Respiratory or skin sensitisation : Not classified

Germ cell mutagenicity : Not classified

Carcinogenicity : Not classified

Reproductive toxicity : Not classified

STOT-single exposure : Not classified

STOT-repeated exposure : Not classified

Aspiration hazard : Not classified

SECTION 12: Ecological information

12.1. Toxicity

Ecology - general : Classification concerning the environment: not applicable.

Ecology - air : Not classified as dangerous for the ozone layer (Regulation (EC) No 1005/2009). TA-Luft Klasse 5.2.1.

Ecology - water : Mild water pollutant (surface water). Maximum concentration in drinking water: 200 mg/l (sodium) (Directive 98/83/EC). Slightly harmful to fishes (LC50(96h) >100 mg/l). Not harmful to invertebrates (Daphnia) (EC50 (48h) > 1000 mg/l). Not harmful to algae (EC50 >1000 mg/l). Not harmful to bacteria (EC50 >1000 mg/l). Nitrification of activated sludge isn't inhibited.

Sodium acetate, trihydrate (6131-90-4)	
LC50 fish 1	100 mg/l (LC50; 964 h)
EC50 Daphnia 1	> 1000 mg/l (EC50; 48 h)

12.2. Persistence and degradability

Sodium acetate, trihydrate (6131-90-4)	
Persistence and degradability	Readily biodegradable in water.

12.3. Bioaccumulative potential

Sodium acetate, trihydrate (6131-90-4)	
BCF fish 1	< 10 (BCF; 72 h; Leuciscus idus)
BCF other aquatic organisms 1	16000 (BCF; 24 h; Chlorella sp.)
BCF other aquatic organisms 2	29100 (BCF; 120 h; Activated sludge)
Bioaccumulative potential	No bioaccumulation data available.

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

No additional information available

Sodium acetate, trihydrate

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 11/12/2016

Doc No: SDS-969.016/1



SECTION 13: Disposal considerations

13.1. Waste treatment methods

Product/Packaging disposal recommendations : Remove waste in accordance with local and/or national regulations. Dissolve or mix with a combustible solvent. Specific preliminary treatment. Remove to an authorized dump. Remove to an authorized incinerator equipped with an afterburner and a flue gas scrubber with energy recovery.

Additional information : LWCA (the Netherlands): KGA category 03. Can be considered as non hazardous waste according to Directive 2008/98/EC.

SECTION 14: Transport information

In accordance with ADR / RID / IMDG / IATA / ADN

ADR	IMDG	IATA	ADN	RID
14.1. UN number				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
14.2. UN proper shipping name				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
14.3. Transport hazard class(es)				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
14.4. Packing group				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
14.5. Environmental hazards				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
No supplementary information available				

14.6. Special precautions for user

- Overland transport

Not applicable

- Transport by sea

Not applicable

- Air transport

Not applicable

- Inland waterway transport

Not applicable

- Rail transport

Not applicable

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

No REACH Annex XVII restrictions

Sodium acetate, trihydrate is not on the REACH Candidate List

Sodium acetate, trihydrate is not on the REACH Annex XIV List

VOC content : 0 %

15.1.2. National regulations

Germany

Sodium acetate, trihydrate

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 11/12/2016

Doc No: SDS-969.016/1



VwVwS Annex reference	: Water hazard class (WGK) 1, low hazard to waters (Classification according to VwVwS, Annex 1 or 2; ID No. 367)
WGK remark	: Classification water polluting in compliance with Verwaltungsvorschrift wassergefährdender Stoffe (VwVwS) of 27 July 2005 (Anhang 2)
12th Ordinance Implementing the Federal Immission Control Act - 12.BImSchV	: Is not subject of the 12. BImSchV (Hazardous Incident Ordinance)

Netherlands

SZW-lijst van kankerverwekkende stoffen	: The substance is not listed
SZW-lijst van mutagene stoffen	: The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding	: The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Vruchtbaarheid	: The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling	: The substance is not listed

15.2. Chemical safety assessment

No additional information available

SECTION 16: Other information

Abbreviations and acronyms:

ADN	European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
ATE	Acute Toxicity Estimate
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
EC50	Median effective concentration
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
PBT	Persistent Bioaccumulative Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
RID	Regulations concerning the International Carriage of Dangerous Goods by Rail
SDS	Safety Data Sheet
vPvB	Very Persistent and Very Bioaccumulative

Data sources : REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

Full text of H- and EUH-statements:

EUH210	Safety data sheet available on request
--------	--

SDS ISOLAB

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product