

SPECIFICATION SHEET

Product Name	:	Sodium Disulfite For analysis ACS,Reag. Ph Eur.
Product Number	:	969.076
Brand	:	ISOLAB
CAS Number	:	7681-57-4
Formula	:	Na ₂ O ₅ S ₂
Formula Weight	:	190.11 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (Na ₂ S ₂ O ₅)	≥ 98.0 - 99.5	%
Identity	passes test	
Appearance of solution	passes test	
pH (10% sol)	4.0 - 5.5	
Sodium Sulfite (Na ₂ SO ₃)	≤ 1.5	
Sodium Sulfate (Na ₂ SO ₄)	≤ 2.0	
Chloride (Cl)	≤ 0.005	%
Heavy metals (as Pb)	≤ 0.001	%
As (Arsenic)	≤ 0.00005	%
Cd (Cadmium)	≤ 0.000005	%
Fe (Iron)	≤ 0.0005	%
Insoluble matter	0	%
Sodium thiosulfate(Na ₂ S ₂ O ₃)	≤ 0.02	%

Quality Management
Isolab Laborgerate GmbH

