



SPECIFICATION SHEET

Product Name : Sodium Disulfite For analysis ACS, Reag. Ph Eur.
Product Number : 969.076
Brand : ISOLAB
CAS Number : 7681-57-4
Formula : $\text{Na}_2\text{O}_5\text{S}_2$
Formula Weight : 190.11 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay ($\text{Na}_2\text{S}_2\text{O}_5$)	≥ 97.0	%
Identity	passes test	
Appearance of solution	passes test	
pH (10% sol)	4.0 - 5.5	
Sodium Sulfite (Na_2SO_3)	≤ 1.5	
Sodium Sulfate (Na_2SO_4)	≤ 2.0	
Chloride (Cl)	≤ 0.005	%
Heavy metals (as Pb)	≤ 0.001	%
As (Arsenic)	≤ 0.00005	%
Cd (Cadmium)	≤ 0.000005	%
Fe (Iron)	≤ 0.0005	%
Insoluble matter	0	%
Sodium thiosulfate ($\text{Na}_2\text{S}_2\text{O}_3$)	≤ 0.02	%

Quality Management
Isolab Laborgerate GmbH

