

# Di-sodium hydrogen phosphate, dodecahydrate

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 05/04/2017 Version: 0.0


Doc. No: SDS-969.096/1



### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

Product form : Substance  
Substance name : di-Sodium hydrogen phosphate dodecahydrate  
EC-No. : 231-448-7  
CAS-No. : 10039-32-4  
Type of product : Pure substance  
Formula : Na<sub>2</sub>HPO<sub>4</sub>·12H<sub>2</sub>O  
Chemical structure :



Synonyms : disodium hydrogen phosphate, 12-hydrate / di-sodium hydrogen phosphate, dodecahydrate / disodium phosphate, 12-hydrate / disodiummonohydrogenphosphate, dodecahydrate / disodiumorthophosphate, dodecahydrate / disodiumphosphate, dodecahydrate / DSP (=disodium phosphate), dodecahydrate / phosphoric acid, disodium salt, 12-hydrate / phosphoric acid, disodium salt, dodecahydrate / secondary-sodium phosphate, dodecahydrate / sec-sodium phosphate, dodecahydrate / sodium hydrogen orthophosphate, dodecahydrate / sodium monohydrogen orthophosphate, dodecahydrate / sodium monohydrogen phosphate, dodecahydrate / sodium orthophosphate dibasic, dodecahydrate / sodium orthophosphate, mono-H, dodecahydrate / sodium phosphate, dibasic, 12-hydrate / sodium phosphate, dibasic, dodecahydrate

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

##### 1.2.1. Relevant identified uses

Use of the substance/mixture : Food industry: additive  
Laboratory chemical

##### 1.2.2. Uses advised against

No additional information available

#### 1.3. Details of the supplier of the safety data sheet

ISOLAB Laborgeräte GmbH  
Am Dillhof 2 - 63863 Eschau / GERMANY  
Tel: + 49 93 74 / 978 55-0  
Fax: +49 93 74 / 978 55-29  
[prodsafe@isolab.de](mailto:prodsafe@isolab.de)

#### 1.4. Emergency telephone number

Country	Organisation/Company	Address	Emergency number	Comment
Germany	Giftnotruf der Charité CBF, Haus VIII (Wirtschaftsgebäude), UG	Hindenburgdamm 30 12203 Berlin	+49 30 19240	

### SECTION 2: Hazards identification

#### 2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not classified

Adverse physicochemical, human health and environmental effects

No additional information available

#### 2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

No labelling applicable

# Di-sodium hydrogen phosphate, dodecahydrate

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 05/04/2017 Version: 0.0

Doc. No: SDS-969.096/1



### 2.3. Other hazards

No additional information available

## SECTION 3: Composition/information on ingredients

### 3.1. Substances

Name	Product identifier	%
di-Sodium hydrogen phosphate dodecahydrate	(CAS-No.) 10039-32-4 (EC-No.) 231-448-7	100

Full text of H-statements: see section 16

### 3.2. Mixtures

Not applicable

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

First-aid measures general	: Check the vital functions. Unconscious: maintain adequate airway and respiration. Respiratory arrest: artificial respiration or oxygen. Cardiac arrest: perform resuscitation. Victim conscious with laboured breathing: half-seated. Victim in shock: on his back with legs slightly raised. Vomiting: prevent asphyxia/aspiration pneumonia. Prevent cooling by covering the victim (no warming up). Keep watching the victim. Give psychological aid. Keep the victim calm, avoid physical strain. Depending on the victim's condition: doctor/hospital.
First-aid measures after inhalation	: Remove the victim into fresh air. Respiratory problems: consult a doctor/medical service.
First-aid measures after skin contact	: Rinse with water. Remove clothing before washing. Take victim to a doctor if irritation persists.
First-aid measures after eye contact	: Rinse with water. Do not apply neutralizing agents. Take victim to an ophthalmologist if irritation persists.
First-aid measures after ingestion	: Rinse mouth with water. Immediately after ingestion: give lots of water to drink. Call Poison Information Centre ( <a href="http://www.big.be/antigif.htm">www.big.be/antigif.htm</a> ). Consult a doctor/medical service if you feel unwell. Ingestion of large quantities: immediately to hospital.

### 4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects after skin contact	: Slight irritation.
Symptoms/effects after eye contact	: Slight irritation.
Symptoms/effects after ingestion	: AFTER ABSORPTION OF HIGH QUANTITIES: Gastrointestinal complaints.
Chronic symptoms	: ON CONTINUOUS/REPEATED EXPOSURE/CONTACT: Red skin. Dry skin.

### 4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

## SECTION 5: Fire-fighting measures

### 5.1. Extinguishing media

Suitable extinguishing media	: Adapt extinguishing media to the environment.
Unsuitable extinguishing media	: No unsuitable extinguishing media known.

### 5.2. Special hazards arising from the substance or mixture

Fire hazard	: DIRECT FIRE HAZARD. No fire hazard. INDIRECT FIRE HAZARD. No fire hazard.
Explosion hazard	: DIRECT EXPLOSION HAZARD. Not applicable. INDIRECT EXPLOSION HAZARD. Not applicable.

### 5.3. Advice for firefighters

Precautionary measures fire	: Exposure to fire/heat: keep upwind. Exposure to fire/heat: consider evacuation. Exposure to fire/heat: have neighbourhood close doors and windows.
Firefighting instructions	: Dilute toxic gases with water spray.
Protection during firefighting	: Heat/fire exposure: compressed air/oxygen apparatus.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

#### 6.1.1. For non-emergency personnel

Protective equipment	: Gloves. Protective clothing. Dust cloud production: compressed air/oxygen apparatus. See "Material-Handling" to select protective clothing.
Emergency procedures	: Mark the danger area. Prevent dust cloud formation, e.g. by wetting. No naked flames. Wash contaminated clothes.

# Di-sodium hydrogen phosphate, dodecahydrate

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 05/04/2017 Version: 0.0

Doc. No: SDS-969.096/1



Measures in case of dust release : In case of dust production: keep upwind. Dust production: have neighbourhood close doors and windows.

### 6.1.2. For emergency responders

No additional information available

### 6.2. Environmental precautions

No additional information available

### 6.3. Methods and material for containment and cleaning up

For containment : Contain released substance, pump into suitable containers. Consult "Material-handling" to select material of containers. Plug the leak, cut off the supply. Knock down/dilute dust cloud with water spray.

Methods for cleaning up : Stop dust cloud by humidifying. Scoop solid spill into closing containers. See "Material-handling" for suitable container materials. Wash clothing and equipment after handling.

### 6.4. Reference to other sections

No additional information available

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

Precautions for safe handling : Comply with the legal requirements. Thoroughly clean/dry the installation before use. Avoid raising dust. Keep away from naked flames/heat. Observe normal hygiene standards. Keep container tightly closed. Carry operations in the open/under local exhaust/ventilation or with respiratory protection.

### 7.2. Conditions for safe storage, including any incompatibilities

Heat and ignition sources : KEEP SUBSTANCE AWAY FROM: heat sources.

Information on mixed storage : KEEP SUBSTANCE AWAY FROM: (strong) acids.

Storage area : Store at ambient temperature. Store in a dry area. Keep container in a well-ventilated place. Meet the legal requirements.

Special rules on packaging : SPECIAL REQUIREMENTS: closing. clean. correctly labelled. meet the legal requirements. Secure fragile packagings in solid containers.

Packaging materials : SUITABLE MATERIAL: paper. textile. polyethylene. cardboard. synthetic material. MATERIAL TO AVOID: aluminium.

### 7.3. Specific end use(s)

No additional information available

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

No additional information available

### 8.2. Exposure controls

#### Materials for protective clothing:

GIVE GOOD RESISTANCE: neoprene. PVC

#### Hand protection:

Gloves

#### Eye protection:

Safety glasses. In case of dust production: protective goggles

#### Skin and body protection:

Protective clothing

#### Respiratory protection:

Dust production: dust mask with filter type P1

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state : Solid

# Di-sodium hydrogen phosphate, dodecahydrate

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 05/04/2017 Version: 0.0

Doc. No: SDS-969.096/1



Appearance	: Crystalline solid. Powder. Grains.
Molecular mass	: 358.14 g/mol
Colour	: Colourless to white.
Odour	: Odourless.
Odour threshold	: No data available
pH	: 9.0 - 9.3 (50 g/l, H <sub>2</sub> O, 20 °C)
pH solution	: 5 %
Relative evaporation rate (butylacetate=1)	: No data available
Melting point	: 35 °C
Freezing point	: No data available
Boiling point	: Not applicable
Flash point	: Not applicable
Auto-ignition temperature	: Not applicable
Decomposition temperature	: > 35 °C
Flammability (solid, gas)	: No data available
Vapour pressure	: No data available
Relative vapour density at 20 °C	: No data available
Relative density	: 1.5
Density	: 1.52 g/cm <sup>3</sup> (20 °C)
Solubility	: Soluble in water. Water: 21.8 g/100ml (20°C)
Log Pow	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive properties	: No data available
Oxidising properties	: No data available
Explosive limits	: No data available

### 9.2. Other information

Minimum ignition energy	: Not applicable
VOC content	: 0 %
Other properties	: Translucent. Substance has basic reaction.

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

On burning: release of toxic and corrosive gases/vapours (phosphorus oxides). Reacts with (some) acids.

### 10.2. Chemical stability

Stable under normal conditions.

### 10.3. Possibility of hazardous reactions

No additional information available

### 10.4. Conditions to avoid

No additional information available

### 10.5. Incompatible materials

No additional information available

### 10.6. Hazardous decomposition products

No additional information available

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

Acute toxicity : Not classified

di-Sodium hydrogen phosphate dodecahydrate (10039-32-4)	
LD50 oral rat	> 2000 mg/kg (Rat)
LD50 dermal rat	> 2000 mg/kg (Rat)

# Di-sodium hydrogen phosphate, dodecahydrate

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 05/04/2017 Version: 0.0

Doc. No: SDS-969.096/1



Skin corrosion/irritation	: Not classified pH: 9 (1 %)
Serious eye damage/irritation	: Not classified pH: 9 (1 %)
Respiratory or skin sensitisation	: Not classified
Germ cell mutagenicity	: Not classified
Carcinogenicity	: Not classified
Reproductive toxicity	: Not classified
STOT-single exposure	: Not classified
STOT-repeated exposure	: Not classified
Aspiration hazard	: Not classified

## SECTION 12: Ecological information

### 12.1. Toxicity

Ecology - air	: Not classified as dangerous for the ozone layer (Regulation (EC) No 1005/2009). TA-Luft Klasse 5.2.1.
Ecology - water	: Mild water pollutant (surface water). Maximum concentration in drinking water: 200 mg/l (sodium) (Directive 98/83/EC). Not harmful to fishes. Slightly harmful to invertebrates (EC50: 100 - 1000 mg/l). May cause eutrophication. Slightly harmful to bacteria. pH shift.

di-Sodium hydrogen phosphate dodecahydrate (10039-32-4)	
LC50 fish 1	> 2400 mg/l (LC50; 48 h)
EC50 Daphnia 1	126 ppm (TLm; 72 h)

### 12.2. Persistence and degradability

di-Sodium hydrogen phosphate dodecahydrate (10039-32-4)	
Persistence and degradability	Biodegradability: not applicable.
Biochemical oxygen demand (BOD)	Not applicable
Chemical oxygen demand (COD)	Not applicable
ThOD	Not applicable

### 12.3. Bioaccumulative potential

di-Sodium hydrogen phosphate dodecahydrate (10039-32-4)	
Bioaccumulative potential	Not bioaccumulative.

### 12.4. Mobility in soil

No additional information available

### 12.5. Results of PBT and vPvB assessment

No additional information available

### 12.6. Other adverse effects

No additional information available

## SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

Product/Packaging disposal recommendations	: Precipitate/make insoluble. Remove to an authorized dump (Class I). Do not discharge into surface water.
Additional information	: Can be considered as non hazardous waste according to Directive 2008/98/EC.
European List of Waste (LoW) code	: 06 03 14 - solid salts and solutions other than those mentioned in 06 03 11 and 06 03 13

## SECTION 14: Transport information

In accordance with ADR / RID / IMDG / IATA / ADN

ADR	IMDG	IATA	ADN	RID
<b>14.1. UN number</b>				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable

# Di-sodium hydrogen phosphate, dodecahydrate

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 05/04/2017 Version: 0.0

Doc. No: SDS-969.096/1



ADR	IMDG	IATA	ADN	RID
<b>14.2. UN proper shipping name</b>				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
<b>14.3. Transport hazard class(es)</b>				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
<b>14.4. Packing group</b>				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
<b>14.5. Environmental hazards</b>				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
No supplementary information available				

### 14.6. Special precautions for user

#### - Overland transport

Not applicable

#### - Transport by sea

Not applicable

#### - Air transport

Not applicable

#### - Inland waterway transport

Not applicable

#### - Rail transport

Not applicable

### 14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### 15.1.1. EU-Regulations

No REACH Annex XVII restrictions

di-Sodium hydrogen phosphate dodecahydrate is not on the REACH Candidate List

di-Sodium hydrogen phosphate dodecahydrate is not on the REACH Annex XIV List

VOC content : 0 %

#### 15.1.2. National regulations

##### Germany

VwVwS Annex reference : Water hazard class (WGK) 1, low hazard to waters (Classification according to VwVwS, Annex 1 or 2; ID No. 330)

WGK remark : Classification water polluting in compliance with Verwaltungsvorschrift wassergefährdender Stoffe (VwVwS) of 27 July 2005 (Anhang 2)

12th Ordinance Implementing the Federal Immission Control Act - 12.BImSchV : Is not subject of the 12. BImSchV (Hazardous Incident Ordinance)

##### Netherlands

SZW-lijst van kankerverwekkende stoffen : The substance is not listed

SZW-lijst van mutagene stoffen : The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding : The substance is not listed

# Di-sodium hydrogen phosphate, dodecahydrate

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 05/04/2017 Version: 0.0

Doc. No: SDS-969.096/1



NIET-limitatieve lijst van voor de voortplanting : The substance is not listed  
giftige stoffen – Vruchtbaarheid

NIET-limitatieve lijst van voor de voortplanting : The substance is not listed  
giftige stoffen – Ontwikkeling

### 15.2. Chemical safety assessment

No additional information available

## SECTION 16: Other information

Abbreviations and acronyms:

ADN	European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
ATE	Acute Toxicity Estimate
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
EC50	Median effective concentration
IARC	International Agency for Research on Cancer
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
PBT	Persistent Bioaccumulative Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
RID	Regulations concerning the International Carriage of Dangerous Goods by Rail
SDS	Safety Data Sheet
vPvB	Very Persistent and Very Bioaccumulative

Data sources : REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

SDS ISOLAB

*This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product*