



SPECIFICATION SHEET

| | | |
|----------------|---|--------------------------------|
| Product Name | : | Sodium hydroxide, pellets pure |
| Product Number | : | 969.112 |
| Brand | : | ISOLAB |
| CAS Number | : | 1310-73-2 |
| Formula | : | NaOH |
| Formula Weight | : | 40.00 g/mol |

| TEST NAME | SPECIFICATION | UNIT |
|---|---------------|------|
| Assay (acidimetric, NaOH) | ≥ 99 % | % |
| Carbonate (as Na ₂ CO ₃) | ≤ 1.0 % | % |
| Chloride (Cl) | ≤ 0.005 % | % |
| Nitrogen Compounds(N) | ≤ 0.0005 % | % |
| Phosphate (PO ₄) | ≤ 0.02 % | % |
| Silicate (SiO ₂) | ≤ 0.01 % | % |
| Sulphate (SO ₄) | ≤ 0.005 % | % |
| Al (Aluminium) | ≤ 0.002 % | % |
| Ca (Calcium) | ≤ 0.005 % | % |
| Copper (Cu) | ≤ 0.0005 % | % |
| Fe (Iron) | ≤ 0.001 % | % |
| Ni (Nickel) | ≤ 0.001 % | % |
| K (Potassium) | ≤ 0.05 % | % |
| Heavy Metals (as Pb) | ≤ 0.0005 % | % |

Quality Management
Isolab Laborgerate GmbH

