



PRODUCT SPECIFICATION

Product Name	:	Sodium Arsenite $\geq 90\%$ for synthesis
Product Number	:	969.271
Brand	:	ISOLAB
CAS Number	:	7784-46-5
Formula	:	AsNaO_2
Formula Weight	:	129.91 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (Iodimetric)	$\geq 90\%$	%
Appearance (Color)	White to yellow and faint Grey to Grey	
Appearance (Form)	Powder	
Solubility (Color)	Colourless to yellow and faint grey to grey and yellow- grey and greyyellow.	
Solubility (Turbidity) 100 mg/mL, H ₂ O	Clear to turbid.	
Antimony	$\leq 0.5\%$	%

Quality Management
Isolab Laborgerate GmbH

