



SPECIFICATION SHEET

Product Name : Starch, Soluble extra pure
Product Number : 970.012
Brand : ISOLAB
CAS Number : 9005-84-9
Formula : $(C_6H_{10}O_5)_n$
Formula Weight : -

TEST NAME	SPECIFICATION	UNIT
pH	6 - 8	
Ash content	≤ 0.5 %	%
Sulphite content (as SO_2)	≤ 0.0005 %	%
Iron (Fe)	≤ 0.001 %	%
Arsenic (As)	≤ 0.000005 %	%
Cadmium (Cd)	≤ 0.000005 %	%
Mercury (Hg)	≤ 0.000005 %	%
Lead (Pb)	≤ 0.00001 %	%
Total aerobic mes. count	≤ 10000 CFU/g	CFU/g
Yeasts	≤ 100 CFU/g	CFU/g
Moulds	≤ 250 CFU/g	CFU/g
Bacillus Cereus	≤ 100 CFU/g	CFU/g
Enterobacteriaceae	≤ 100 CFU/g	CFU/g
Escherichia coli in 1 g	absent	
Salmonella in 25 g	absent	

Quality Management
Isolab Laborgerate GmbH

