



PRODUCT SPECIFICATION

Product Name : Sulfuric Acid 95-98% , Extra pure
Product Number : 970.022
Brand : ISOLAB
CAS Number : 7664-93-9
Formula : H_2SO_4
Formula Weight : 98.08 g/mol

| TEST NAME | SPECIFICATION | UNIT |
|------------|----------------------------|----------|
| Assay | 95 – 98 | % |
| Density | ~ 1,84 | g/cm^3 |
| Solubility | Soluble in water | |
| Colour | Colorless, Slightly blurry | |

Quality Management
Isolab Laborgerate GmbH

