

# TARTRAZINE

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 21/02/2017


Doc No: SDS-971.D01/2



### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

Product form	: Substance
Substance name	: TARTRAZINE AR FOR MICROSCOPY
EC No	: 217-699-5
CAS No	: 1934-21-0
Type of product	: Pure substance
Formula	: C <sub>16</sub> H <sub>9</sub> N <sub>4</sub> Na <sub>3</sub> O <sub>9</sub> S <sub>2</sub>
Chemical structure	:



Synonyms	: 1310 yellow / 1409 yellow / 1H-pyrazole-3-carboxylic acid, 4,5-dihydro-5-oxo-1-(4-sulfophenyl)-4-[(4-sulfophenyl)azo]-, trisodium salt / 3-carboxy-5-hydroxy-1-p-sulfophenyl-4-p-sulfophenylazopyrazole trisodium salt / 4,5-dihydro-5-oxo-1-(4-sulfophenyl)-4-[(4-sulfophenyl)azo]-1H-pyrazole-3-carboxylic acid trisodium salt / 5-hydroxy-1-(para-sulfophenyl)-4-[(para-sulfophenyl)azo]pyrazole-3-carboxylic acid trisodium salt / 5-hydroxy-1-(p-sulfophenyl)-4-(p-sulfophenyl)azopyrazole-3-carboxylic acid, trisodium salt / 5-hydroxy-1-(p-sulfophenyl)-4-[(p-sulfophenyl)azo]pyrazole-3-carboxylic acid trisodium salt / A.F. yellow No. 4 / acid leather yellow T / acid yellow 23 / acid yellow T / acilan yellow GG / airedale yellow T / aizen tartrazine / amacid yellow T / atul tartrazine / bucadid tartrazine / C.I. 19140 / C.I. 640 / C.I. acid yellow 23 / C.I. acid yellow 23, trisodium salt / C.I. food yellow 4 / calcocid yellow MCG / calcocid yellow XX / canacert tartrazine / certicol tartrazol yellow S / cilefa yellow T / curon fast yellow 5G / D and C yellow No. 5 / dolkwal tartrazine / dye FD and C yellow No. 5 / E 102 / edicol supra tartrazine N / egg yellow A / erio tartrazine / eurocert tartrazine / FD & C yellow No. 5 / FD & C yellow No. 5 tartrazine / FD and C yellow 5 / fenazo yellow T / food yellow 4 / food yellow 5 / food yellow No. 4 / HD tartrazine / HD tartrazine supra / hexacert yellow No. 5 / hexacol tartrazine / hidazid tartrazine / hispacid fast yellow T / hydrazine yellow / hydroxine yellow L / kako tartrazine / kayaku food colour yellow No. 4 / kayaku tartrazine / KCA foodcol tartrazine PF / KCA tartrazine PF / kiton yellow T / lake yellow / lemon yellow A / lemon yellow A GEIGY / maple tartrazol yellow / mitsui tartrazine / naphtocard yellow O / neklacid yellow T / oxanal yellow T / pyrazole-3-carboxylic acid, 5-hydroxy-1-(p-sulfophenyl)-4-(p-sulfophenyl)azo-, trisodium salt / san-ei tartrazine / Schultz No. 737 / sugai tartrazine / tartar yellow FS / tartar yellow N / tartar yellow PF / tartar yellow S / tartran yellow / tartraphenine / tartrazine A expo T / tartrazine B / tartrazine B.P.C. / tartrazine extra pure A / tartrazine FD & C yellow #5 / tartrazine FQ / tartrazine G / tartrazine lake / tartrazine lake yellow N / tartrazine M / tartrazine MCGL / tartrazine N / tartrazine NS / tartrazine O / tartrazine T / tartrazine XX / tartrazine XXX / tartrazine yellow / tartrazol BPC / tartrazol yellow / tartrine yellow O / trisodium 3-carboxy-5-hydroxy-1-p-sulfophenyl-4-p-sulfophenylazopyrazole / trisodium 5-hydroxy-1-(4-sulfonatophenyl)-4-(4-sulfonatophenylazo)pyrazole-3-carboxylate / trisodium 5-hydroxy-1-(4-sulphophenyl)-4-(4-sulphophenylazo)pyrazole-3-carboxylate / trisodium salt of 3-carboxy-5-hydroxy-1-sulfophenylazopyrazole / unitertracid yellow TE / usacert yellow No. 5 / vondacid tartrazine / wool yellow / xylene fast yellow GT / Y-4 / yellow lake 69 / yellow No. 5 / yellow No. 5 FDC
----------	---

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

##### 1.2.1. Relevant identified uses

Industrial/Professional use spec	: For professional use only
Use of the substance/mixture	: Food industry: dyestuff Cosmetic product: dyestuff Laboratory chemical

##### 1.2.2. Uses advised against

No additional information available

#### 1.3. Details of the supplier of the safety data sheet

ISOLAB GmbH

# TARTRAZINE

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 2/21/2017

Doc No: SDS-971.D01/2



ISOLAB Laborgeräte GmbH  
Am Dillhof 2 - 63863 Eschau / GERMANY  
Tel: + 49 93 74 / 978 55-0  
Fax: +49 93 74 / 978 55-29  
[prodsafe@isolab.de](mailto:prodsafe@isolab.de)

### 1.4. Emergency telephone number

Country	Organisation/Company	Address	Emergency number	Comment
Germany	Giftnotruf der Charité CBF, Haus VIII (Wirtschaftgebäude), UG	Hindenburgdamm 30 12203 Berlin	+49 30 19240	

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

#### Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin. Sens. 1 H317  
Eye Irrit. 2 H319  
Resp. Sens. 1 H334

Full text of hazard classes and H-statements : see section 16

#### Adverse physicochemical, human health and environmental effects

No additional information available

### 2.2. Label elements

#### Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP)



Signal word (CLP)

: Danger

Hazard statements (CLP)

: H319 - Causes serious eye irritation.  
H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled.  
H317 - May cause an allergic skin reaction.

Precautionary statements (CLP)

: P285 - In case of inadequate ventilation wear respiratory protection.  
P280 - Wear protective gloves, protective clothing, eye protection, face protection.  
P261 - Avoid breathing dust, fume, gas, mist, vapours, spray.  
P264 - Wash thoroughly after handling.  
P272 - Contaminated work clothing should not be allowed out of the workplace  
P302+P352 - IF ON SKIN : Wash with plenty of soap and water.  
P305+P351+P338 - IF IN EYES : Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P304+P341 - IF INHALED : If breathing is difficult, remove to fresh air and keep at rest in a position comfortable for breathing.  
P333+P313 - If skin irritation or rash occurs : Get medical advice.  
P342+P311 - If experiencing respiratory symptoms : Call a POISON CENTER or doctor.  
P363 - Wash contaminated clothing before reuse.  
P501 - Dispose of this material and its container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation.

### 2.3. Other hazards

This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII

This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII

## SECTION 3: Composition/information on ingredients

### 3.1. Substances

Name	Product identifier	%
TARTRAZINE	(CAS No) 1934-21-0 (EC No) 217-699-5	100

Full text of H-statements: see section 16

# TARTRAZINE

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 2/21/2017

Doc No: SDS-971.D01/2



### 3.2. Mixtures

Not applicable

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

First-aid measures general	: Check the vital functions. Unconscious: maintain adequate airway and respiration. Respiratory arrest: artificial respiration or oxygen. Cardiac arrest: perform resuscitation. Victim conscious with laboured breathing: half-seated. Victim in shock: on his back with legs slightly raised. Vomiting: prevent asphyxia/aspiration pneumonia. Prevent cooling by covering the victim (no warming up). Keep watching the victim. Give psychological aid. Keep the victim calm, avoid physical strain. Depending on the victim's condition: doctor/hospital.
First-aid measures after inhalation	: Remove the victim into fresh air. Respiratory problems: consult a doctor/medical service.
First-aid measures after skin contact	: Rinse with water. Soap may be used. Take victim to a doctor if irritation persists.
First-aid measures after eye contact	: Rinse with water. Do not apply neutralizing agents. Take victim to an ophthalmologist if irritation persists.
First-aid measures after ingestion	: Rinse mouth with water. Give activated charcoal. Call Poison Information Centre ( <a href="http://www.big.be/antigif.htm">www.big.be/antigif.htm</a> ). Consult a doctor/medical service if you feel unwell. Ingestion of large quantities: immediately to hospital.

### 4.2. Most important symptoms and effects, both acute and delayed

Chronic symptoms	: Causes serious eye irritation. May cause allergy or asthma symptoms or breathing difficulties if inhaled. May cause an allergic skin reaction.
------------------	--

### 4.3. Indication of any immediate medical attention and special treatment needed

General information	: Never give anything by mouth to an unconscious person. If you feel unwell, seek medical advice (show the label where possible).
---------------------	---

## SECTION 5: Fire-fighting measures

### 5.1. Extinguishing media

Suitable extinguishing media	: Water spray. Polyvalent foam. ABC powder. Carbon dioxide.
Unsuitable extinguishing media	: No unsuitable extinguishing media known.

### 5.2. Special hazards arising from the substance or mixture

Fire hazard	: DIRECT FIRE HAZARD. No data available on direct fire hazard. INDIRECT FIRE HAZARD. No data available on indirect fire hazard.
Explosion hazard	: DIRECT EXPLOSION HAZARD. No data available on direct explosion hazard. INDIRECT EXPLOSION HAZARD. No data available on indirect explosion hazard.

### 5.3. Advice for firefighters

Precautionary measures fire	: Exposure to fire/heat: keep upwind. Exposure to fire/heat: consider evacuation. Exposure to fire/heat: have neighbourhood close doors and windows.
Firefighting instructions	: Dilute toxic gases with water spray.
Protection during firefighting	: Heat/fire exposure: compressed air/oxygen apparatus.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

#### 6.1.1. For non-emergency personnel

Protective equipment	: Gloves. Protective clothing. Dust cloud production: compressed air/oxygen apparatus. See "Material-Handling" to select protective clothing.
Emergency procedures	: Mark the danger area. Prevent dust cloud formation, e.g. by wetting. No naked flames. Wash contaminated clothes.

#### 6.1.2. For emergency responders

No additional information available

### 6.2. Environmental precautions

No additional information available

### 6.3. Methods and material for containment and cleaning up

For containment	: Contain released substance, pump into suitable containers. Plug the leak, cut off the supply. Knock down/dilute dust cloud with water spray.
-----------------	--

# TARTRAZINE

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 2/21/2017

Doc No: SDS-971.D01/2



Methods for cleaning up : Prevent dust cloud formation. Scoop solid spill into closing containers. Clean contaminated surfaces with an excess of water. Wash clothing and equipment after handling.

### 6.4. Reference to other sections

No additional information available

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

Precautions for safe handling : Comply with the legal requirements. Remove contaminated clothing immediately. Clean contaminated clothing. Thoroughly clean/dry the installation before use. Avoid raising dust. Keep away from naked flames/heat. Reduce/avoid exposure and/or contact. Keep container tightly closed. Carry operations in the open/under local exhaust/ventilation or with respiratory protection.

### 7.2. Conditions for safe storage, including any incompatibilities

Heat and ignition sources : KEEP SUBSTANCE AWAY FROM: heat sources.  
Information on mixed storage : KEEP SUBSTANCE AWAY FROM: oxidizing agents.  
Storage area : Store in a dry area. Meet the legal requirements.  
Special rules on packaging : SPECIAL REQUIREMENTS: closing. watertight. dry. correctly labelled. meet the legal requirements. Secure fragile packagings in solid containers.  
Packaging materials : SUITABLE MATERIAL: No data available. MATERIAL TO AVOID: No data available.

### 7.3. Specific end use(s)

No additional information available

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

No additional information available

### 8.2. Exposure controls

#### Personal protective equipment:

Gloves. Safety glasses. Dust formation: dust mask.

#### Materials for protective clothing:

GIVE GOOD RESISTANCE: rubber

#### Hand protection:

Gloves

#### Eye protection:

Safety glasses. In case of dust production: protective goggles

#### Skin and body protection:

Protective clothing

#### Respiratory protection:

Dust formation: dust mask



## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state : Solid  
Appearance : Solid. Powder.  
Molecular mass : 534.36 g/mol  
Colour : Orange.  
Odour : Odourless.  
Odour threshold : No data available

# TARTRAZINE

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 2/21/2017

Doc No: SDS-971.D01/2



pH	: No data available
Relative evaporation rate (butylacetate=1)	: No data available
Melting point	: No data available
Freezing point	: No data available
Boiling point	: No data available
Flash point	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Flammability (solid, gas)	: No data available
Vapour pressure	: No data available
Relative vapour density at 20 °C	: No data available
Relative density	: No data available
Solubility	: No data available
Log Pow	: -10.17 (Estimated value)
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive properties	: No data available
Oxidising properties	: No data available
Explosive limits	: No data available

### 9.2. Other information

VOC content	: 0 %
Other properties	: Hygroscopic.

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

On heating/burning: release of toxic and corrosive gases/vapours (nitrous vapours, sulphur oxides, carbon monoxide - carbon dioxide).

### 10.2. Chemical stability

Hygroscopic.

### 10.3. Possibility of hazardous reactions

No additional information available

### 10.4. Conditions to avoid

No additional information available

### 10.5. Incompatible materials

No additional information available

### 10.6. Hazardous decomposition products

No additional information available

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

Acute toxicity	: Not classified
Skin corrosion/irritation	: Not classified
Serious eye damage/irritation	: Not classified
Respiratory or skin sensitisation	: Not classified
Germ cell mutagenicity	: Not classified
Carcinogenicity	: Not classified
Reproductive toxicity	: Not classified
STOT-single exposure	: Not classified
STOT-repeated exposure	: Not classified
Aspiration hazard	: Not classified

# TARTRAZINE

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 2/21/2017

Doc No: SDS-971.D01/2



### SECTION 12: Ecological information

#### 12.1. Toxicity

Ecology - air : Not classified as dangerous for the ozone layer (Regulation (EC) No 1005/2009).  
Ecology - water : Mild water pollutant (surface water). Maximum concentration in drinking water: 200 mg/l (sodium) (Directive 98/83/EC). Not harmful to fishes (LC50 >1000 mg/l).

TARTRAZINE (1934-21-0)	
LC50 fish 1	> 1000 mg/l (LC50; 48 h)

#### 12.2. Persistence and degradability

TARTRAZINE (1934-21-0)	
Persistence and degradability	Biodegradability in water: no data available. Biodegradability in soil: no data available.

#### 12.3. Bioaccumulative potential

TARTRAZINE (1934-21-0)	
BCF fish 1	< 3 (BCF)
Log Pow	-10.17 (Estimated value)
Bioaccumulative potential	Low potential for bioaccumulation (BCF < 500).

#### 12.4. Mobility in soil

No additional information available

#### 12.5. Results of PBT and vPvB assessment

TARTRAZINE (1934-21-0)	
This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII	
This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII	

#### 12.6. Other adverse effects

No additional information available

### SECTION 13: Disposal considerations

#### 13.1. Waste treatment methods

Product/Packaging disposal recommendations : Dissolve or mix with a combustible solvent. Remove to an authorized incinerator equipped with an afterburner and a flue gas scrubber with energy recovery.  
Additional information : Can be considered as non hazardous waste according to Directive 2008/98/EC.  
European List of Waste (LoW) code : 07 03 99 - wastes not otherwise specified

### SECTION 14: Transport information

In accordance with ADR / RID / IMDG / IATA / ADN

ADR	IMDG	IATA	ADN	RID
<b>14.1. UN number</b>				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
<b>14.2. UN proper shipping name</b>				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
<b>14.3. Transport hazard class(es)</b>				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
<b>14.4. Packing group</b>				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
<b>14.5. Environmental hazards</b>				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
No supplementary information available				

#### 14.6. Special precautions for user

##### - Overland transport

Not applicable

# TARTRAZINE

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 2/21/2017

Doc No: SDS-971.D01/2



### - Transport by sea

Not applicable

### - Air transport

Not applicable

### - Inland waterway transport

Not applicable

### - Rail transport

Not applicable

### 14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### 15.1.1. EU-Regulations

No REACH Annex XVII restrictions

TARTRAZINE is not on the REACH Candidate List

TARTRAZINE is not on the REACH Annex XIV List

VOC content : 0 %

#### 15.1.2. National regulations

##### Germany

VwVwS Annex reference : Water hazard class (WGK) 1, low hazard to waters (Classification according to VwVwS, Annex 3; ID No. 3249)

WGK remark : Classification water polluting based on the R-phrases in compliance with Verwaltungsvorschrift wassergefährdender Stoffe (VwVwS) of 27 July 2005 (Anhang 3)

12th Ordinance Implementing the Federal Immission Control Act - 12.BImSchV : Is not subject of the 12. BImSchV (Hazardous Incident Ordinance)

##### Netherlands

SZW-lijst van kankerverwekkende stoffen : The substance is not listed

SZW-lijst van mutagene stoffen : The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding : The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Vruchtbaarheid : The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling : The substance is not listed

### 15.2. Chemical safety assessment

No additional information available

## SECTION 16: Other information

Abbreviations and acronyms:

ADN	European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
ATE	Acute Toxicity Estimate
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
EC50	Median effective concentration
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
PBT	Persistent Bioaccumulative Toxic

# TARTRAZINE

## Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 2/21/2017

Doc No: SDS-971.D01/2



REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
RID	Regulations concerning the International Carriage of Dangerous Goods by Rail
SDS	Safety Data Sheet
vPvB	Very Persistent and Very Bioaccumulative
Skin. Sens. 1	Skin sensitisation - Category 1
Eye Irrit. 2	Eye irritation - Category 2
Resp. Sens. 1	Respiratory sensitization - Category 1
H319	Causes serious eye irritation.
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H317	May cause an allergic skin reaction.

Data sources : REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

SDS ISOLAB

*This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product*