



**SPECIFICATION SHEET**

**Product Name** : Tetrabutylammonium hydroxide 40% in methanol  
**Product Number** : 972.06P  
**Brand** : ISOLAB  
**CAS Number** : -  
**Formula** : -  
**Formula Weight** : -

TEST NAME	SPECIFICATION	UNIT
Description	Clear to slightly hazy, colourless to slight yellow liquid with characteristic odour	
Solubility	Miscible with water	
Assay	≥ 40.00 %	%

Quality Management  
Isolab Laborgerate GmbH

