



SPECIFICATION SHEET

Product Name : Tetraethyl orthosilicate, for synthesis
Product Number : 972.091
Brand : ISOLAB
CAS Number : 78-10-4
Formula : $\text{Si}(\text{OC}_2\text{H}_5)_4$
Formula Weight : 208.33 g/mol

TEST NAME	SPECIFICATION	UNIT
Appearance (Color)	Colorless	
Appearance (Form)	Liquid	
Infrared spectrum	Conforms to Structure	
Purity (GC)	≥ 97.5 %	%

Quality Management
Isolab Laborgerate GmbH

