



SPECIFICATION SHEET

Product Name : Toluene, for liquid chromatography
Product Number : 973.017
Brand : ISOLAB
CAS Number : 108-88-3
Formula : C₇H₈
Formula Weight : 92.14 g/mol

TEST NAME	SPECIFICATION	UNIT
Purity (GC)	≥ 99.9 %	%
Identity (IR)	passes test	
Evaporation residue	≤ 0.0002 %	%
Water	≤ 0.05 %	%
Acidity	≤ 0.0002 meq/g	meq/g
Alkalinity	≤ 0.0006 meq/g	meq/g
Transmission (at 300 nm)	≥ 70 %	%
Transmission (at 310 nm)	≥ 80 %	%
Transmission (from 350 nm)	≥ 98 %	%
Filtered by 0.2 µm filter		

Quality Management
Isolab Laborgerate GmbH

