



SPECIFICATION SHEET

Product Name : Triethylamine, for synthesis
Product Number : 974.031
Brand : ISOLAB
CAS Number : 121-44-8
Formula : $C_6H_{15}N$
Formula Weight : 101.19 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (GC, area%)	$\geq 99.0 \%$ (a/a)	%
Identity (IR)	passes test	
Water (K.F.)	$\leq 0.2 \%$	%

Quality Management
Isolab Laborgerate GmbH

