



SPECIFICATION SHEET

Product Name : Water, for liquid chromatography
Product Number : 983.017
Brand : ISOLAB
CAS Number : 7732-18-5
Formula : H₂O
Formula Weight : 18.02 g/mol

TEST NAME	SPECIFICATION	UNIT
Non-volatile matter	≤ 0.0003 %	%
Fluorescence (chinin) at 254 nm	≤ 0.0000001 %	%
Fluorescence (chinin) at 365 nm	≤ 0.0000001 %	%
HPLC Gradient at 254 nm	≤ 1 mAU	mAU
HPLC Gradient at 210 nm	≤ 5 mAU	mAU
Chloride (Cl)	≤ 0.000001 %	%
Fluoride (F)	≤ 0.000001 %	%
Nitrate (NO ₃)	≤ 0.00001 %	%
Sulfate (SO ₄)	≤ 0.00001 %	%
Filtered via 0.2 µm pore size filter		

Quality Management
Isolab Laborgerate GmbH

