



SPECIFICATION SHEET

Product Name : Ninhydrin, GR for analysis ACS, Reag. Ph Eur
Product Number : 951.036
Brand : ISOLAB
CAS Number : 485-47-2
Formula : $C_9H_6O_4$
Formula Weight : 178.15 g/mol

TEST NAME	SPECIFICATION	UNIT
Identity (IR-spectrum)	passes test	
Appearance	White to yellow or brownish white	%
Solubility (10 g/l, water)	passes test	
Melting point	about 255 °C	°C
Suitability for amino acid analysis	passes test	

Quality Management
Isolab Laborgerate GmbH

